HASELNUSS

Molekulare Allergiediagnostik

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Literatur: 1. Flinterman AE et al. *Clinical reactivity to hazelnut in children: association with sensitization to birch pollen or nuts?* J Allergy Clin Immunol. 2006 Nov; 118(5): 1186–9. 2. De Knop KJ et al. *Age-related sensitization profiles for hazelnut (Corylus avellana) in a birch-endemic region.* Pediatr Allergy Immunol. 2011 Feb; 22(1 Pt 2): e139–49. 3. Hartz C et al. *Comparison of IgE-Binding Capacity, Cross-Reactivity and Biological Potency of Allergenic Non-Specific Lipid Transfer Proteins from Peach, Cherry and Hazelnut.* Int Arch Allergy Immunol. 2010 Jun 17; 153(4): 335–46. 4. Schulten V et al. *Pru p 3, the nonspecific lipid transfer protein from peach, dominates the immune response to its homolog in hazelnut.* Allergy. 2011 66 1005–13. 5. Beyer K et al. *Identification of an 11S globulin as a major hazelnut food allergen in hazelnut-induced systemic reactions.* J Allergy Clin Immunol. 2002 Sep; 110(3): 517–23. 6. Garino C et al. *Isolation, cloning, and characterization of the 2S albumin: a new allergen from hazelnut.* Mol Nutr Food Res. 2010 Sep; 54(9): 1257–65. 7. Masthoff L et al. *Sensitization to Cor a 9 and Cor a 14 is highly specific for a severe hazelnut allergy in Dutch children and adults.* J Allergy Clin Immunol. 2013(In press). 8. Verweij MM et al. *Young infants with atopic dermatitis can display sensitization to Cor a 9, an 11S legumin-like seed-storage protein from hazelnut (Corylus avellana).* Pediatric Allergy and Immunology. 2011; 22(2): 196–201. 9. Hansen KS et al. *Component-resolved in vitro diagnosis of hazelnut allergy in Europe.* J Allergy Clin Immunol. 2009 Apr 1; 123(5): 1134–41. 10. Pastorello EA et al. *Identification of hazelnut major allergens in sensitive patients with positive double-blind, placebo-controlled food challenge results.* J Allergy Clin Immunol. 2002; 109(3): 563–70. 11. Hansen KS et al. *Roasted hazelnuts-allergenic activity evaluated by double-blind, placebo-controlled food challenge.* Allergy. 2003 Feb; 58(2): 132–8. 12. Asero R et al. *Walnut-induced anaphylaxis with cross-reactivity to hazelnut and Brazil nut.* J Allergy Clin Immunol. 2004 Feb; 113(2): 358–60. 13. Flinterman AE et al. *Hazelnut allergy: from pollen-associated mild allergy to severe anaphylactic reactions.* Curr Opin Allergy Clin Immunol. 2008 Jun; 8(3): 261–5. 14. Sicherer SH et al. *US prevalence of self-reported peanut, tree nut, and sesame allergy: 11-year follow-up.* J Allergy Clin Immunol. 2010 May 10; 125(5): 1077–83. 15. Osterballe M et al. *The prevalence of food hypersensitivity in an unselected population of children and adults.* Pediatr Allergy Immunol. 2005 Nov; 16(7): 567–73. 16. Roehr CC et al. *Food allergy and non-allergic food hypersensitivity in children and adolescents.* Clin Exp Allergy. 2004 Oct; 34(10): 1534–41. 17. Schocker F et al. *Recombinant lipid transfer protein Cor a 8 from hazelnut: a new tool for in vitro diagnosis of potentially severe hazelnut allergy.* J Allergy Clin Immunol. 2004 Jan; 113(1): 141–7.

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© 2013 Thermo Fisher Scientific Inc. Alle Rechte vorbehalten. Alle Warenzeichen sind das Eigentum von Thermo Fisher Scientific Inc. und seiner Tochtergesellschaften. Rechtmäßiger Hersteller: Phadia AB, Uppsala, Schweden
Phadia GmbH, Munzinger Str. 7, 79111 Freiburg / Deutschland, Tel. +49 761 478050
Phadia Austria GmbH, Donau-City-Str. 1, 1220 Wien / Österreich, Tel. +43 1 2702020
Phadia AG, Sennweidstr. 46, 6312 Steinhausen / Schweiz, Tel. +41 43 3434050