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# Literatur

1. Valenta R, Lidholm J, Niederberger V, Hayek B, Kraft D, Gronlund H. **The recombinant allergen-based concept of component-resolved diagnostics and immunotherapy (CRD and CRIT).** *Clin Exp Allergy.* 1999 Jul;29(7):896-904.
2. Canonica GW, Ansotegui IJ, Pawankar R, Schmid-Grendelmeier P, van Hage M, Baena-Cagnani CE, et al. **A WAO - ARIA - GA(2)LEN consensus document on molecular-based allergy diagnostics.** *World Allergy Organ J.* 2013;6(1):17.
3. Hiller R, Laffer S, Harwanegg C, Huber M, Schmidt WM, Twardosz A, et al. **Microarrayed allergen molecules: diagnostic gatekeepers for allergy treatment.** *FASEB J.* 2002 Mar;16(3):414-6.
4. Harwanegg C, Hiller R. **Protein microarrays for the diagnosis of allergic diseases: state-of-the-art and future development.** *Clin Chem Lab Med.* 2005;43(12):1321-6.
5. Ferrer M, Sanz ML, Sastre J, Bartra J, del Cuvillo A, Montoro J, et al. **Molecular diagnosis in allergology: application of the microarray technique.** *J Investig Allergol Clin Immunol.* 2009;19 Suppl 1:19-24.
6. Cabrera-Freitag P, Goikoetxea MJ, Gamboa PM, Martinez-Aranguren R, Beorlegui C, Fernandez J, et al. **A study of the variability of the in vitro component-based microarray ISAC CDR 103 technique.** *J Investig Allergol Clin Immunol.* 2011;21(5):414-5.
7. Gadisseur R, Chapelle JP, Cavalier E. **A new tool in the field of in-vitro diagnosis of allergy: preliminary results in the comparison of ImmunoCAP 250 with the ImmunoCAP ISAC.** *Clin Chem Lab Med.* 2011 Feb;49(2):277-80.
8. Lupinek C, Wollmann E, Baar A, Banerjee S, Breiteneder H, Broecker BM, et al. **Advances in allergen-microarray technology for diagnosis and monitoring of allergy: the MeDALL allergen-chip.** *Methods.* 2014 Mar 1;66(1):106-19.
9. Westman M, Lupinek C, Bousquet J, Andersson N, Pahr S, Baar A, et al. **Early childhood IgE reactivity to pathogenesis-related class 10 proteins predicts allergic rhinitis in adolescence.** *J Allergy Clin Immunol.* 2015 May;135(5):1199-206 e1-11.
10. Skrinko I, Lupinek C, Valenta R, Hovland V, Pahr S, Baar A, et al. **The use of the MeDALL-chip to assess IgE sensitization: a new diagnostic tool for allergic disease?** *Pediatr Allergy Immunol.* 2015 May;26(3):239-46.
11. Martinez-Aranguren R, Lizaso MT, Goikoetxea MJ, Garcia BE, Cabrera-Freitag P, Trellez O, et al. **Is the determination of specific IgE against components using ISAC 112 a reproducible technique?** *PLoS One.* 2014;9(2):e88394.
12. Choi JS, Roh JY, Lee JR. **Clinical availability of component-resolved diagnosis using microarray technology in atopic dermatitis.** *Ann Dermatol.* 2014 Aug;26(4):437-46.
13. Lizaso MT, Garcia BE, Tabar AI, Lasa E, Echechipia S, Alvarez MJ, et al. **Comparison of conventional and component-resolved diagnostics by two different methods (Advia-Centaur/Microarray-ISAC) in pollen allergy.** *Ann Allergy Asthma Immunol.* 2011 Jul;107(1):35-41.
14. Huss-Marp PDJ, Gutermuth J, Schäffner I, Darsow U, Pfab F, Brockow K, et al. **Comparison of molecular and extract-based allergy diagnostics with multiplex and singleplex analysis.** *Allergo Journal International.* 2015;24(2):46-53.
15. Melioli G, Bonifazi F, Bonini S, Maggi E, Mussap M, Passalacqua G, et al. **The ImmunoCAP ISAC molecular allergology approach in adult multi-sensitized Italian patients with respiratory symptoms.** *Clin Biochem.* 2011 Aug;44(12):1005-11.
16. Scala E, Alessandri C, Bernardi ML, Ferrara R, Palazzo P, Pomponi D, et al. **Cross-sectional survey on immunoglobulin E reactivity in 23,077 subjects using an allergenic molecule-based microarray detection system.** *Clin Exp Allergy.* 2010 Jun;40(6):911-21.
17. Rossi RE, Melioli G, Monasterolo G, Harwanegg C, Rossi L, Canonica GW, et al. **Sensitization profiles in polysensitized patients from a restricted geographical area: further lessons from multiplexed component resolved diagnosis.** *Eur Ann Allergy Clin Immunol.* 2011 Dec;43(6):171-5.
18. Burney PG, Potts J, Kummeling I, Mills EN, Clausen M, Dubakiene R, et al. **The prevalence and distribution of food sensitization in European adults.** *Allergy.* 2014 Mar;69(3):365-71.
19. Panzner P, Vachova M, Vitovcova P, Brodska P, Vlas T. **A comprehensive analysis of middle-European molecular sensitization profiles to pollen allergens.** *Int Arch Allergy Immunol.* 2014;164(1):74-82.
20. Passalacqua G, Melioli G, Bonifazi F, Bonini S, Maggi E, Senna G, et al. **The additional values of micro-array allergen assay in the management of polysensitized patients with respiratory allergy.** *Allergy.* 2013 Aug;68(8):1029-33.
21. Alvarado MI, Jimeno L, De La Torre F, Boissy P, Rivas B, Lazaro MJ, et al. **Profilin as a severe food allergen in allergic patients overexposed to grass pollen.** *Allergy.* 2014 Dec;69(12):1610-6.

22. van Kampen V, Sander I, Quirce S, Bruning T, Merget R, Raulf M. **IgE Sensitization to Lupine in Bakers - Cross-Reactivity or Co-Sensitization to Wheat Flour?** *Int Arch Allergy Immunol.* 2015;166(1):63-70.
23. Ebo DG, Bridts CH, Verweij MM, De Knop KJ, Hagedorens MM, De Clerck LS, et al. **Sensitization profiles in birch pollen-allergic patients with and without oral allergy syndrome to apple: lessons from multiplexed component-resolved allergy diagnosis.** *Clin Exp Allergy.* 2010 Feb;40(2):339-47.
24. Patelis A, Gunnbjornsdottir M, Malinovsky A, Matsson P, Onell A, Hogman M, et al. **Population-based study of multiplexed IgE sensitization in relation to asthma, exhaled nitric oxide, and bronchial responsiveness.** *J Allergy Clin Immunol.* 2012 Aug;130(2):397-402 e2.
25. Nordlund B, Konradsen JR, Kull I, Borres MP, Onell A, Hedlin G, et al. **IgE antibodies to animal-derived lipocalin, kallikrein and secretoglobulin are markers of bronchial inflammation in severe childhood asthma.** *Allergy.* 2012 May;67(5):661-9.
26. Konradsen JR, Nordlund B, Onell A, Borres MP, Gronlund H, Hedlin G. **Severe childhood asthma and allergy to furry animals: refined assessment using molecular-based allergy diagnostics.** *Pediatr Allergy Immunol.* 2014 Mar;25(2):187-92.
27. Bjerg A, Winberg A, Berthold M, Mattsson L, Borres MP, Ronmark E. **A population-based study of animal component sensitization, asthma, and rhinitis in schoolchildren.** *Pediatr Allergy Immunol.* 2015 Jun 9. DOI: 10.1111/pai.12422.
28. Gray CL, Levin ME, Zar HJ, Potter PC, Khumalo NP, Volkwyn L, et al. **Food allergy in South African children with atopic dermatitis.** *Pediatr Allergy Immunol.* 2014 Oct;25(6):572-9.
29. Fung I, Kim JS, Spergel JM. **Relating microarray component testing and reported food allergy and food-triggered atopic dermatitis: a real-world analysis.** *Ann Allergy Asthma Immunol.* 2013 Mar;110(3):173-7 e1.
30. Ott H, Weissmantel S, Kennes LN, Merk HF, Baron JM, Folster-Holst R. **Molecular microarray analysis reveals allergen- and exotoxin-specific IgE repertoires in children with atopic dermatitis.** *J Eur Acad Dermatol Venereol.* 2014 Jan;28(1):100-7.
31. Rockmann H, van Geel MJ, Knulst AC, Huiskes J, Bruijn-zeel-Koomen CA, de Bruin-Weller MS. **Food allergen sensitization pattern in adults in relation to severity of atopic dermatitis.** *Clin Transl Allergy.* 2014;4(1):9.
32. Mari A, Scala E, Alessandri C. **The IgE-microarray testing in atopic dermatitis: a suitable modern tool for the immunological and clinical phenotyping of the disease.** *Curr Opin Allergy Clin Immunol.* 2011 Oct;11(5):438-44.
33. Assa'ad AH, Putnam PE, Collins MH, Akers RM, Jameson SC, Kirby CL, et al. **Pediatric patients with eosinophilic esophagitis: an 8-year follow-up.** *J Allergy Clin Immunol.* 2007 Mar;119(3):731-8.
34. Plaza-Martin AM, Jimenez-Feijoo R, Andaluz C, Giner-Munoz MT, Martin-Mateos MA, Piquer-Gibert M, et al. **Polysensitization to aeroallergens and food in eosinophilic esophagitis in a pediatric population.** *Allergol Immunopathol (Madr).* 2007 Jan-Feb;35(1):35-7.
35. Armentia A, Martin S, Barrio J, Martin B, Garcia JC, Vega JM, et al. **Value of microarray allergen assay in the management of eosinophilic oesophagitis.** *Allergol Immunopathol (Madr).* 2015 Jan-Feb;43(1):73-80.
36. Erwin EA, James HR, Gutekunst HM, Russo JM, Kelleher KJ, Platts-Mills TA. **Serum IgE measurement and detection of food allergy in pediatric patients with eosinophilic esophagitis.** *Ann Allergy Asthma Immunol.* 2010 Jun;104(6):496-502.
37. Simon D, Straumann A, Dahinden C, Simon HU. **Frequent sensitization to *Candida albicans* and profilins in adult eosinophilic esophagitis.** *Allergy.* 2013 Jul;68(7):945-8.
38. van Rhijn BD, van Ree R, Versteeg SA, Vlieg-Boerstra BJ, Sprikkelman AB, Terreehorst I, et al. **Birch pollen sensitization with cross-reactivity to food allergens predominates in adults with eosinophilic esophagitis.** *Allergy.* 2013 Nov;68(11):1475-81.
39. Ackerbauer D, Bublun M, Radauer C, Varga EM, Hafner C, Ebner C, et al. **Component-Resolved IgE Profiles in Austrian Patients with a Convincing History of Peanut Allergy.** *Int Arch Allergy Immunol.* 2015 Feb 27;166(1):13-24.
40. Klemans RJ, Liu X, Knulst AC, Knol MJ, Gmelig-Meyling F, Borst E, et al. **IgE binding to peanut components by four different techniques: Ara h 2 is the most relevant in peanut allergic children and adults.** *Clin Exp Allergy.* 2013 Aug;43(8):967-74.
41. D'Urbano LE, Pellegrino K, Artesani MC, Donnanno S, Luciano R, Riccardi C, et al. **Performance of a component-based allergen-microarray in the diagnosis of cow's milk and hen's egg allergy.** *Clin Exp Allergy.* 2010 Oct;40(10):1561-70.
42. Ebo DG, Hagedorens MM, De Knop KJ, Verweij MM, Bridts CH, De Clerck LS, et al. **Component-resolved diagnosis from latex allergy by microarray.** *Clin Exp Allergy.* 2010 Feb;40(2):348-58.
43. Seyfarth F, Schliemann S, Wiegand C, Hipler UC, Elsner P. **Diagnostic value of the ISAC allergy chip in detecting latex sensitizations.** *Int Arch Occup Environ Health.* 2014 Oct;87(7):775-81.

44. Hong X, Caruso D, Kumar R, Liu R, Liu X, Wang G, et al. **IgE, but not IgG4, antibodies to Ara h 2 distinguish peanut allergy from asymptomatic peanut sensitization.** *Allergy*. 2012 Dec;67(12):1538-46.
45. Ciprandi G, Pistorio A, Silvestri M, Rossi GA, Tosca MA. **Peanut anaphylaxis: the usefulness of molecular-based allergy diagnostics.** *Allergy*. 2015 Jan;70(1):129-30.
46. Schuler S, Ferrari G, Schmid-Grendelmeier P, Harr T. **Microarray-based component-resolved diagnosis of latex allergy: isolated IgE-mediated sensitization to latexprofilin Hev b8 may act as confounder.** *Clin Transl Allergy*. 2013;3(1):11.
47. Gassner M, Gehrig R, Schmid-Grendelmeier P. **Hay Fever as a Christmas Gift.** *N Engl J Med*. 2013;368(4):393-4.
48. Cardona V, Luengo O, Garriga T, Labrador-Horrillo M, Sala-Cunill A, Izquierdo A, et al. **Co-factor-enhanced food allergy.** *Allergy*. 2012 Oct;67(10):1316-8.
49. Pascal M, Munoz-Cano R, Reina Z, Palacin A, Vilella R, Picado C, et al. **Lipid transfer protein syndrome: clinical pattern, cofactor effect and profile of molecular sensitization to plant-foods and pollens.** *Clin Exp Allergy*. 2012 Oct;42(10):1529-39.
50. Romano A, Scala E, Rumi G, Gaeta F, Caruso C, Alonzi C, et al. **Lipid transfer proteins: the most frequent sensitizer in Italian subjects with food-dependent exercise-induced anaphylaxis.** *Clin Exp Allergy*. 2012 Nov;42(11):1643-53.
51. Carter S, Heaps A, Boswijk K, Jolles S, Kaminski E. **Identification of clinically relevant allergens using the Phadia ISAC microarray in patients with idiopathic anaphylaxis.** *Clin Exp Allergy*. 2012;42(12):1829-30.
52. Heaps A, Carter S, Selwood C, Moody M, Unsworth J, Deacock S, et al. **The utility of the ISAC allergen array in the investigation of idiopathic anaphylaxis.** *Clin Exp Immunol*. 2014 Aug;177(2):483-90.
53. Letran A, Espinazo M, Moreno F. **Measurement of IgE to pollen allergen components is helpful in selecting patients for immunotherapy.** *Ann Allergy Asthma Immunol*. 2013 Oct;111(4):295-7.
54. Sastre J, Landivar ME, Ruiz-Garcia M, Andregnette-Rosigno MV, Mahillo I. **How molecular diagnosis can change allergen-specific immunotherapy prescription in a complex pollen area.** *Allergy*. 2012 May;67(5):709-11.
55. Asero R. **Component-resolved diagnosis-assisted prescription of allergen-specific immunotherapy: a practical guide.** *Eur Ann Allergy Clin Immunol*. 2012 Oct;44(5):183-7.
56. Wollmann E, Lupinek C, Kundi M, Selb R, Niederberger V, Valenta R. **Reduction in allergen-specific IgE binding as measured by microarray: A possible surrogate marker for effects of specific immunotherapy.** *J Allergy Clin Immunol*. 2015 Apr 23. DOI: 10.1016/j.jaci.2015.02.034.
57. Hatzler L, Panetta V, Lau S, Wagner P, Bergmann RL, Illi S, et al. **Molecular spreading and predictive value of preclinical IgE response to Phleum pratense in children with hay fever.** *J Allergy Clin Immunol*. 2012 Oct;130(4):894-901 e5.
58. Onell A, Hjalte L, Borres MP. **Exploring the temporal development of childhood IgE profiles to allergen components.** *Clin Transl Allergy*. 2012;2(1):24.
59. Prosperi MC, Belgrave D, Buchan I, Simpson A, Custovic A. **Challenges in interpreting allergen microarrays in relation to clinical symptoms: a machine learning approach.** *Pediatr Allergy Immunol*. 2014 Feb;25(1):71-9.
60. Simpson A, Lazic N, Belgrave DC, Johnson P, Bishop C, Mills C, et al. **Patterns of IgE responses to multiple allergen components and clinical symptoms at age 11 years.** *J Allergy Clin Immunol*. 2015 Apr 30. DOI: 10.1016/j.jaci.2015.03.027.
61. Melioli G, Marcomini L, Agazzi A, Bazurro G, Tosca M, Rossi GA, et al. **The IgE repertoire in children and adolescents resolved at component level: a cross-sectional study.** *Pediatr Allergy Immunol*. 2012 Aug;23(5):433-40.

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© 2017 Thermo Fisher Scientific Inc. Alle Rechte vorbehalten. Alle Warenzeichen sind das Eigentum von Thermo Fisher Scientific und seiner Tochtergesellschaften, falls nicht anders angegeben. Rechtmäßiger Hersteller: Phadia AB, Uppsala, Schweden
**Phadia GmbH**, Munzinger Str. 7, D-79111 Freiburg, Tel. +49 761 47 8050, Fax +49 761 47 805 338
**Thermo Fisher Diagnostics Austria GmbH**, Dresdner Str. 89, A-1200 Wien, Tel. +43 1 270 20 20, Fax +43 1 270 20 20 20
**Phadia AG**, Sennweidstr. 46, CH-6312 Steinhausen, Tel. +41 43 343 40 50, Fax +41 43 343 40 51