**MOLEKULARE  
ALLERGIEDIAGNOSTIK**

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1. Zuidmeer L, van Ree R. Lipid transfer protein allergy: primary food allergy or pollen/food syndrome in some cases. *Curr Opin Allergy Clin Immunol.* 2007;7:269-73.
2. Richard C, Leduc V, Battais F. Plant lipid transfer proteins (LTPs): biochemical aspect in panallergen – structural and functional features, and allergenicity. *Allerg Immunol (Paris).* 2007;39:76-84.
3. Wopfner N, Dissertori O, Ferreira F, Lackner P. Calcium-binding proteins and their role in allergic diseases. *Immunol Allergy Clin North Am.* 2007;27:29-44.
4. Bohle B. The impact of pollen-related food allergens on pollen allergy. *Allergy.* 2007;62:3-10.
5. Salcedo G, Sanchez-Monge R, Barber D, Diaz-Perales A. Plant non-specific lipid transfer proteins: an interface between plant defence and human allergy. *Biochim Biophys Acta.* 2007;1771:781-91.
6. Marion D, Bakan B, Elmorjani K. Plant lipid binding proteins: properties and applications. *Biotechnol Adv.* 2007;25:195-7.
7. Altmann F. The role of protein glycosylation in allergy. *Int Arch Allergy Immunol.* 2007;142:99-115.
8. Egger M, Mutschlechner S, Wopfner N, Gadermaier G, Briza P, Ferreira F. Pollen-food syndromes associated with weed pollinosis: an update from the molecular point of view. *Allergy.* 2006;61:461-76.
9. Breiteneder H, Clare Mills EN. Plant food allergens – structural and functional aspects of allergenicity. *Biotechnol Adv.* 2005;23:395-9.
10. Gangur V, Kelly C, Navuluri L. Sesame allergy: a growing food allergy of global proportions? *Ann Allergy Asthma Immunol.* 2005;95:4-11.
11. Breiteneder H, Mills C. Nonspecific lipid-transfer proteins in plant foods and pollens: an important allergen class. *Curr Opin Allergy Clin Immunol.* 2005;5:275-9.
12. Breiteneder H, Mills EN. Molecular properties of food allergens. *J Allergy Clin Immunol.* 2005;115:14-23.
13. Malandain H. IgE-reactive carbohydrate epitopes – classification, cross-reactivity, and clinical impact. *Allerg Immunol (Paris).* 2005;37:122-8.
14. Salcedo G, Sanchez-Monge R, Diaz-Perales A, Garcia-Casado G, Barber D. Plant non-specific lipid transfer proteins as food and pollen allergens. *Clin Exp Allergy.* 2004;34:1336-41.
15. Mothes N, Horak F, Valenta R. Transition from a botanical to a molecular classification in tree pollen allergy: implications for diagnosis and therapy. *Int Arch Allergy Immunol.* 2004;135:357-73.
16. Mills EN, Jenkins JA, Alcocer MJ, Shewry PR. Structural, biological, and evolutionary relationships of plant food allergens sensitizing via the gastrointestinal tract. *Crit Rev Food Sci Nutr.* 2004;44:379-407.
17. Breiteneder H, Radauer C. A classification of plant food allergens. *J Allergy Clin Immunol.* 2004;113:821-30.
18. Roux KH, Teuber SS, Sathe SK. Tree nut allergens. *Int Arch Allergy Immunol.* 2003;131:234-44.
19. Wagner S, Breiteneder H. The latex-fruit syndrome. *Biochem Soc Trans.* 2002;30:935-40.
20. Hoffmann-Sommergruber K. Pathogenesis-related (PR)-proteins identified as allergens. *Biochem Soc Trans.* 2002;30:930-5.
21. Shewry PR, Beaudoin F, Jenkins J, Griffiths-Jones S, Mills EN. Plant protein families and their relationships to food allergy. *Biochem Soc Trans.* 2002 ;30:906-10.
22. van Ree R. Clinical importance of non-specific lipid transfer proteins as food allergens. *Biochem Soc Trans.* 2002;30:910-3.
23. Yagami T. Allergies to cross-reactive plant proteins. Latex-fruit syndrome is comparable with pollen-food allergy syndrome. *Int Arch Allergy Immunol.* 2002;128:271-9.
24. Vieths S, Scheurer S, Ballmer-Weber B. Current understanding of cross-reactivity of food allergens and pollen. *Ann N Y Acad Sci.* 2002;964:47-68.
25. Midoro-Horiuti T, Brooks EG, Goldblum RM. Pathogenesis-related proteins of plants as allergens. *Ann Allergy Asthma Immunol.* 2001;87:261-71.
26. Sicherer SH. Clinical implications of cross-reactive food allergens. *J Allergy Clin Immunol.* 2001;108:881-90.
27. Aalberse RC, Akkerdaas J, van Ree R. Cross-reactivity of IgE antibodies to allergens. *Allergy.* 2001;56:478-90.
28. Pastorello EA, Pompei C, Pravettoni V, Brenna O, Farioli L, Trambaioli C, Conti A. Lipid transfer proteins and 2S albumins as allergens. *Allergy.* 2001;56 Suppl 67:45-7. Review.
29. Salcedo G, Diaz-Perales A, Sanchez-Monge R. The role of plant panallergens in sensitization to natural rubber latex. *Curr Opin Allergy Clin Immunol.* 2001;1:177-83.
30. Breiteneder H, Ebner C. Molecular and biochemical classification of plant-derived food allergens. *J Allergy Clin Immunol.* 2000;106:27-36.
31. Caballero T, Martin-Esteban M. Association between pollen hypersensitivity and edible vegetable allergy: a review. *J Investig Allergol Clin Immunol.* 1998;8:6-16.

## Literatur zu PR-10 Proteinen (Bet v 1-Homologen)

1. Ballmer-Weber BK, Holzhauser T, Scibilia J, Mittag D, Zisa G, Ortolani C, Oesterballe M, Poulsen LK, Vieths S, Bindslev-Jensen C. Clinical characteristics of soybean allergy in Europe: a double-blind, placebo-controlled food challenge study. *J Allergy Clin Immunol.* 2007;119:1489-96.
2. Bohle B. The impact of pollen-related food allergens on pollen allergy. *Allergy.* 2007;62:3-10. Review.
3. Moverare R, Kosunen TU, Haahtela T. Change in the pattern of IgE reactivity to timothy grass and birch pollen allergens over a 20-year period. *J Investig Allergol Clin Immunol.* 2006;16:274-8.
4. Fernandez-Rivas M, Bolhaar S, Gonzalez-Mancebo E, Asero R, van Leeuwen A, Bohle B, Ma Y, Ebner C, Rigby N, Sancho AI, Miles S, Zuidmeer L, Knulst A, Breiteneder H, Mills C, Hoffmann-Sommergruber K, van Ree R. Apple allergy across Europe: how allergen sensitization profiles determine the clinical expression of allergies to plant foods. *J Allergy Clin Immunol.* 2006;118:481-8.
5. Bohle B, Zwolfer B, Heratizadeh A, Jahn-Schmid B, Antonia YD, Alter M, Keller W, Zuidmeer L, van Ree R, Werfel T, Ebner C. Cooking birch pollen-related food: divergent consequences for IgE- and T cell-mediated reactivity in vitro and in vivo. *J Allergy Clin Immunol.* 2006;118:242-9.

6. Asero R, Marzban G, Martinelli A, Zaccarini M, Machado ML. Search for low-allergenic apple cultivars for birch-pollen-allergic patients: is there a correlation between in vitro assays and patient response? *Allerg Immunol (Paris)*. 2006;38:94-8.
7. Moreno-Ancillo A, Gil-Adrados AC, Cosmes PM, Dominguez-Noche C, Pineda F. Role of Dau c 1 in three different patterns of carrot-induced asthma. *Allergol Immunopathol (Madr)*. 2006;34:116-20.
8. Cudowska B, Kaczmarek M, Restani P. Immunoblotting in the diagnosis of cross-reactivity in children allergic to birch. *Rocz Akad Med Bialymst*. 2005;50:268-73.
9. Schimek EM, Zwolfer B, Briza P, Jahn-Schmid B, Vogel L, Vieths S, Ebner C, Bohle B. Gastrointestinal digestion of Bet v 1-homologous food allergens destroys their mediator-releasing, but not T cell-activating, capacity. *J Allergy Clin Immunol*. 2005;116:1327-33.
10. Bohle B, Radakovics A, Luttkopf D, Jahn-Schmid B, Vieths S, Ebner C. Characterization of the T cell response to the major hazelnut allergen, Cor a 1.04: evidence for a relevant T cell epitope not cross-reactive with homologous pollen allergens. *Clin Exp Allergy*. 2005;35:1392-9.
11. Mittag D, Batori V, Neudecker P, Wiche R, Friis EP, Ballmer-Weber BK, Vieths S, Roggen EL. A novel approach for investigation of specific and cross-reactive IgE epitopes on Bet v 1 and homologous food allergens in individual patients. *Mol Immunol*. 2006;43:268-78.
12. Mittag D, Vieths S, Vogel L, Wagner-Loew D, Starke A, Hunziker P, Becker WM, Ballmer-Weber BK. Birch pollen-related food allergy to legumes: identification and characterization of the Bet v 1 homologue in mungbean (*Vigna radiata*), Vig r 1. *Clin Exp Allergy*. 2005;35:1049-55.
13. Schirmer T, Hoffmann-Sommergrube K, Susani M, Breiteneder H, Markovic-Housley Z. Crystal structure of the major celery allergen Api g 1: molecular analysis of cross-reactivity. *J Mol Biol*. 2005;351:1101-9.
14. Ballmer-Weber BK, Wangorsch A, Bohle B, Kaul S, Kundig T, Fotisch K, van Ree R, Vieths S. Component-resolved in vitro diagnosis in carrot allergy: does the use of recombinant carrot allergens improve the reliability of the diagnostic procedure? *Clin Exp Allergy*. 2005;35:970-8.
15. Jahn-Schmid B, Radakovics A, Luttkopf D, Scheurer S, Vieths S, Ebner C, Bohle B. Bet v 1142-156 is the dominant T-cell epitope of the major birch pollen allergen and important for cross-reactivity with Bet v 1-related food allergens. *J Allergy Clin Immunol*. 2005;116:213-9.
16. Ricci G, Righetti F, Menna G, Bellini F, Miniaci A, Masi M. Relationship between Bet v 1 and Bet v 2 specific IgE and food allergy in children with grass pollen respiratory allergy. *Mol Immunol*. 2005;42:1251-7.
17. Purohit A, Laffer S, Metz-Favre C, Verot A, Kricek F, Valenta R, Pauli G. Poor association between allergen-specific serum immunoglobulin E levels, skin sensitivity and basophil degranulation: a study with recombinant birch pollen allergen Bet v 1 and an immunoglobulin E detection system measuring immunoglobulin E capable of binding to Fc epsilon RI. *Clin Exp Allergy*. 2005;35:186-92.
18. Erdmann SM, Sachs B, Schmidt A, Merk HF, Scheiner O, Moll-Stodowy S, Sauer I, Kwiecien R, Maderegger B, Hoffmann-Sommergruber K. In vitro analysis of birch-pollen-associated food allergy by use of recombinant allergens in the basophil activation test. *Int Arch Allergy Immunol*. 2005;136:230-8.
19. Ghunaim N, Gronlund H, Kronqvist M, Gronneberg R, Soderstrom L, Ahlstedt S, van Hage-Hamsten M. Antibody profiles and self-reported symptoms to pollen-related food allergens in grass pollen-allergic patients from northern Europe. *Allergy*. 2005;60:185-91.
20. Cudowska B, Kaczmarek M. Diagnostic value of birch recombinant allergens (rBet v 1, profilin rBet v 2) in children with pollen-related food allergy. *Rocz Akad Med Bialymst*. 2004;49:111-5.
21. Bolhaar ST, van Ree R, Ma Y, Bruijnzeel-Koomen CA, Vieths S, Hoffmann-Sommergruber K, Knulst AC, Zuidmeer L. Severe allergy to sharon fruit caused by birch pollen. *Int Arch Allergy Immunol*. 2005;136:45-52.
22. Mittag D, Akkerdaas J, Ballmer-Weber BK, Vogel L, Wensing M, Becker WM, Koppelman SJ, Knulst AC, Helbling A, Hefle SL, Van Ree R, Vieths S. Ara h 8, a Bet v 1-homologous allergen from peanut, is a major allergen in patients with combined birch pollen and peanut allergy. *J Allergy Clin Immunol*. 2004;114:1410-7.
23. Karlsson AL, Alm R, Ekstrand B, Fjellner-Modig S, Schiott A, Bengtsson U, Bjork L, Hjerno K, Roepstorff P, Emanuelsson CS. Bet v 1 homologues in strawberry identified as IgE-binding proteins and presumptive allergens. *Allergy*. 2004;59:1277-84.
24. Scheurer S, Lauer I, Foetisch K, San Miguel Moncin M, Retzek M, Hartz C, Enrique E, Lidholm J, Cistero-Bahima A, Vieths S. Strong allergenicity of Pru av 3, the lipid transfer protein from cherry, is related to high stability against thermal processing and digestion. *J Allergy Clin Immunol*. 2004;114:900-7.
25. Bolhaar ST, Ree R, Bruijnzeel-Koomen CA, Knulst AC, Zuidmeer L. Allergy to jackfruit: a novel example of Bet v 1-related food allergy. *Allergy*. 2004;59:1187-92.
26. Breuer K, Wulf A, Constien A, Tetau D, Kapp A, Werfel T. Birch pollen-related food as a provocation factor of allergic symptoms in children with atopic eczema/dermatitis syndrome. *Allergy*. 2004;59:988-94.
27. Mittag D, Vieths S, Vogel L, Becker WM, Rihs HP, Helbling A, Wutrich B, Ballmer-Weber BK. Soybean allergy in patients allergic to birch pollen: clinical investigation and molecular characterization of allergens. *J Allergy Clin Immunol*. 2004;113:148-54.
28. De Amici M, Mosca M, Vignini M, Quagliani S, Moratti R. Recombinant birch allergens (Bet v 1 and Bet v 2) and the oral allergy syndrome in patients allergic to birch pollen. *Ann Allergy Asthma Immunol*. 2003;91:490-2.
29. Bohle B, Radakovics A, Jahn-Schmid B, Hoffmann-Sommergruber K, Fischer GF, Ebner C. Bet v 1, the major birch pollen allergen, initiates sensitization to Api g 1, the major allergen in celery: evidence at the T cell level. *Eur J Immunol*. 2003;33:3303-10.
30. Fernandez-Rivas M, Gonzalez-Mancebo E, Rodriguez-Perez R, Benito C, Sanchez-Monge R, Salcedo G, Alonso MD, Rosado A, Tejedor MA, Vila C, Casas ML. Clinically relevant peach allergy is related to peach lipid transfer protein, Pru p 3, in the Spanish population. *J Allergy Clin Immunol*. 2003;112:789-95.
31. Mari A, Wallner M, Ferreira F. Fagales pollen sensitization in a birch-free area: a respiratory cohort survey using Fagales pollen extracts and birch recombinant allergens (rBet v 1, rBet v 2, rBet v 4). *Clin Exp Allergy*. 2003;33:1419-28.
32. Asero R, Mistrello G, Roncarolo D, Amato S, Zanoni D, Barocci F, Caldironi G. Detection of clinical markers of sensitization to profilin in patients allergic to plant-derived foods. *J Allergy Clin Immunol*. 2003;112:427-32.

33. Cadot P, Kochuyt AM, van Ree R, Ceuppens JL. Oral allergy syndrome to chicory associated with birch pollen allergy. *Int Arch Allergy Immunol.* 2003;131:19-24.
34. Bodtger U, Poulsen LK, Malling HJ. Asymptomatic skin sensitization to birch predicts later development of birch pollen allergy in adults: a 3-year follow-up study. *J Allergy Clin Immunol.* 2003;111:149-54.
35. Moverare R, Westritschnig K, Svensson M, Hayek B, Bende M, Pauli G, Sorva R, Haahtela T, Valenta R, Elfman L. Different IgE reactivity profiles in birch pollen-sensitive patients from six European populations revealed by recombinant allergens: an imprint of local sensitization. *Int Arch Allergy Immunol.* 2002;128:325-35.
36. Wensing M, Akkerdaas JH, van Leeuwen WA, Stapel SO, Bruijnzeel-Koomen CA, Aalberse RC, Bast BJ, Knulst AC, van Ree R. IgE to Bet v 1 and profilin: cross-reactivity patterns and clinical relevance. *J Allergy Clin Immunol.* 2002;110:435-42.
37. Vieths S, Luttkopf D, Reindl J, Anliker MD, Wuthrich B, Ballmer-Weber BK. Allergens in celery and zucchini. *Allergy.* 2002;57 Suppl 72:100-5.
38. Ballmer-Weber BK, Scheurer S, Fritsche P, Enrique E, Cistero-Bahima A, Haase T, Wuthrich B. Component-resolved diagnosis with recombinant allergens in patients with cherry allergy. *J Allergy Clin Immunol.* 2002;110:167-73.
39. Pastorello EA, Vieths S, Pravettoni V, Farioli L, Trambaioli C, Fortunato D, Luttkopf D, Calamari M, Ansaloni R, Scibilia J, Ballmer-Weber BK, Poulsen LK, Wutrich B, Hansen KS, Robino AM, Ortolani C, Conti A. Identification of hazelnut major allergens in sensitive patients with positive double-blind, placebo-controlled food challenge results. *J Allergy Clin Immunol.* 2002;109:563-70.
40. Ballmer-Weber BK, Wuthrich B, Wangorsch A, Fotisch K, Altmann F, Vieths S. Carrot allergy: double-blinded, placebo-controlled food challenge and identification of allergens. *J Allergy Clin Immunol.* 2001;108:301-7.
41. Anhoj C, Backer V, Nolte H. Diagnostic evaluation of grass- and birch-allergic patients with oral allergy syndrome. *Allergy.* 2001;56:548-52.
42. Holm J, Baerentzen G, Gajhede M, Ipsen H, Larsen JN, Lowenstein H, Wissenbach M, Spangfort MD. Molecular basis of allergic cross-reactivity between group 1 major allergens from birch and apple. *J Chromatogr B Biomed Sci Appl.* 2001;756:307-13.
43. Karamloo F, Scheurer S, Wangorsch A, May S, Hausteiner D, Vieths S. Pyr c 1, the major allergen from pear (*Pyrus communis*), is a new member of the Bet v 1 allergen family. *J Chromatogr B Biomed Sci Appl.* 2001;756:281-93.
44. Mari A. Multiple pollen sensitization: a molecular approach to the diagnosis. *Int Arch Allergy Immunol.* 2001;125:57-65.
45. Paschke A, Kinder H, Zunker K, Wigotzki M, Steinhart H, Wessbecher R, Vieluf I. Characterization of cross-reacting allergens in mango fruit. *Allergy.* 2001;56:237-42.
46. Reider N, Sepp N, Fritsch P, Weinlich G, Jensen-Jarolim E. Anaphylaxis to camomile: clinical features and allergen cross-reactivity. *Clin Exp Allergy.* 2000;30:1436-43.
47. Hoffmann-Sommergruber K, Demoly P, Cramer R, Breiteneder H, Ebner C, Laimer Da Camara Machado M, Blaser K, Ismail C, Scheiner O, Bousquet J, Menz G. IgE reactivity to Api g 1, a major celery allergen, in a Central European population is based on primary sensitization by Bet v 1. *J Allergy Clin Immunol.* 1999;104:478-84.
48. Reekers R, Busche M, Wittmann M, Kapp A, Werfel T. Birch pollen-related foods trigger atopic dermatitis in patients with specific cutaneous T-cell responses to birch pollen antigens. *J Allergy Clin Immunol.* 1999;104:466-72.
49. Moneo I, Gomez M, Sanchez-Monge R, Alday E, de las Heras M, Esteban I, Bootello A, Salcedo G. Lack of crossreaction with Bet v 1 in patients sensitized to Dau c 1, a carrot allergen. *Ann Allergy Asthma Immunol.* 1999;83:71-5.
50. Werfel T, Reekers R, Busche M, Schmidt P, Constien A, Wittmann M, Kapp A. Association of birch pollen-related food-responsive atopic dermatitis with birch pollen allergen-specific T-cell reactions. *Curr Probl Dermatol.* 1999;28:18-28.
51. Van Ree R, Van Leeuwen WA, Akkerdaas JH, Aalberse RC. How far can we simplify in vitro diagnostics for Fagales tree pollen allergy? A study with three whole pollen extracts and purified natural and recombinant allergens. *Clin Exp Allergy.* 1999;29:848-55.
52. Hoffmann-Sommergruber K, O'Riordain G, Ahorn H, Ebner C, Laimer Da Camara Machado M, Puhlinger H, Scheiner O, Breiteneder H. Molecular characterization of Dau c 1, the Bet v 1 homologous protein from carrot and its cross-reactivity with Bet v 1 and Api g 1. *Clin Exp Allergy.* 1999;29:840-7.
53. Ebner C, Jensen-Jarolim E, Leitner A, Breiteneder H. Characterization of allergens in plant-derived spices: Apiaceae spices, pepper (Piperaceae), and paprika (bell peppers, Solanaceae). *Allergy.* 1998;53(46 Suppl):52-4.
54. Fritsch R, Bohle B, Vollmann U, Wiedermann U, Jahn-Schmid B, Krebitz M, Breiteneder H, Kraft D, Ebner C. Bet v 1, the major birch pollen allergen, and Mal d 1, the major apple allergen, cross-react at the level of allergen-specific T helper cells. *J Allergy Clin Immunol.* 1998;102:679-86.
55. Vega A, Dominguez C, Cosmes P, Martinez A, Bartolome B, Martinez J, Palacios R. Anaphylactic reaction to ingestion of *Quercus ilex* acorn nut. *Clin Exp Allergy.* 1998;28:739-42.
56. Jensen-Jarolim E, Santner B, Leitner A, Grimm R, Scheiner O, Ebner C, Breiteneder H. Bell peppers (*Capsicum annuum*) express allergens (profilin, pathogenesis-related protein P23 and Bet v 1) depending on the horticultural strain. *Int Arch Allergy Immunol.* 1998;116:103-9.
57. Leitner A, Jensen-Jarolim E, Grimm R, Wuthrich B, Ebner H, Scheiner O, Kraft D, Ebner C. Allergens in pepper and paprika. Immunologic investigation of the celery-birch-mugwort-spice syndrome. *Allergy.* 1998;53:36-41.
58. Fernandez-Rivas M, van Ree R, Cuevas M. Allergy to Rosaceae fruits without related pollinosis. *J Allergy Clin Immunol.* 1997;100:728-33.
59. Jensen-Jarolim E, Leitner A, Hirschwehr R, Kraft D, Wuthrich B, Scheiner O, Graf J, Ebner C. Characterization of allergens in Apiaceae spices: anise, fennel, coriander and cumin. *Clin Exp Allergy.* 1997;27:1299-306.
60. Elfman L, Svensson M, Lidholm J, Pauli G, Valenta R. Different profiles in specific IgE to rBet v 1 and rBet v 2 in patients allergic to birch pollen from six countries. *Int Arch Allergy Immunol.* 1997;113:249-51.
61. Rossi RE, Monasterolo G, Operti D, Corsi M. Evaluation of recombinant allergens Bet v 1 and Bet v 2 (profilin) by Pharmacia CAP system in patients with pollen-related allergy to birch and apple. *Allergy.* 1996;51:940-5.

62. Pauli G, Oster JP, Deviller P, Heiss S, Bessot JC, Susani M, Ferreira F, Kraft D, Valenta R. Skin testing with recombinant allergens rBet v 1 and birch profilin, rBet v 2: diagnostic value for birch pollen and associated allergies. *J Allergy Clin Immunol.* 1996;97:1100-9.
63. Ebner C, Hirschehr R, Bauer L, Breiteneder H, Valenta R, Ebner H, Kraft D, Scheiner O. Identification of allergens in fruits and vegetables: IgE cross-reactivities with the important birch pollen allergens Bet v 1 and Bet v 2 (birch profilin). *J Allergy Clin Immunol.* 1995;95:962-9.
64. Ebner C, Birkner T, Valenta R, Rumpold H, Breitenbach M, Scheiner O, Kraft D. Common epitopes of birch pollen and apples – studies by western and northern blot. *J Allergy Clin Immunol.* 1991;88:588-94.
11. Vassilopoulou E, Zuidmeer L, Akkerdaas J, Tassios I, Rigby NR, Mills EN, van Ree R, Saxoni-Papageorgiou P, Papadopoulos NG. Severe immediate allergic reactions to grapes: part of a lipid transfer protein-associated clinical syndrome. *Int Arch Allergy Immunol.* 2007;143:92-102.
12. Asero R, Amato S, Alfieri B, Folloni S, Mistrello G. Rice: another potential cause of food allergy in patients sensitized to lipid transfer protein. *Int Arch Allergy Immunol.* 2007;143:69-74.
13. Primavesi L, Brenna OV, Pompei C, Pravettoni V, Farioli L, Pastorello EA. Influence of cultivar and processing on cherry (*Prunus avium*) allergenicity. *J Agric Food Chem.* 2006;54:9930-5.
14. Rihs HP, Rueff F, Lundberg M, Rozynek P, Barber D, Scheurer S, Cistero-Bahima A, Bruning T, Raulf-Heimsoth M. Relevance of the recombinant lipid transfer protein of *Hevea brasiliensis*: IgE-binding reactivity in fruit-allergic adults. *Ann Allergy Asthma Immunol.* 2006;97:643-9.

## Literatur zu nsLTP (PR-14 Proteinen)

1. Cudowska B, Kaczmarek M, Restani P. Lipid transfer protein in diagnosis of birch-apple syndrome in children. *Immunobiology.* 2008;213:89-96.
2. Flinterman AE, Akkerdaas JH, den Hartog Jager CF, Rigby NM, Fernandez-Rivas M, Hoekstra MO, Bruijnzeel-Koomen CA, Knulst AC, van Ree R, Pasmans SG. Lipid transfer protein-linked hazelnut allergy in children from a non-Mediterranean birch-endemic area. *J Allergy Clin Immunol.* 2008 Feb;121:423-28.
3. Borges JP, Barre A, Culerrier R, Granier C, Didier A, Rougé P. Lipid transfer proteins from Rosaceae fruits share consensus epitopes responsible for their IgE-binding cross-reactivity. *Biochem Biophys Res Commun.* 2008 Jan 25;365:685-90.
4. Palacin A, Quirce S, Armentia A, Fernández-Nieto M, Pacios LF, Asensio T, Sastre J, Diaz-Perales A, Salcedo G. Wheat lipid transfer protein is a major allergen associated with baker's asthma. *J Allergy Clin Immunol.* 2007 Nov;120:1132-8.
5. Asero R, Mistrello G, Roncarolo D, Amato S. Detection of some safe plant-derived foods for LTP-allergic patients. *Int Arch Allergy Immunol.* 2007;144:57-63.
6. Pastorello EA, Farioli L, Conti A, Pravettoni V, Bonomi S, Iametti S, Fortunato D, Scibilia J, Bindslev-Jensen C, Ballmer-Weber B, Robino AM, Ortolani C. Wheat IgE-mediated food allergy in European patients: alpha-amylase inhibitors, lipid transfer proteins and low-molecular-weight glutenins. Allergenic molecules recognized by double-blind, placebo-controlled food challenge. *Int Arch Allergy Immunol.* 2007;144:10-22.
7. Gamboa PM, Caceres O, Antepara I, Sanchez-Monge R, Ahrazem O, Salcedo G, Barber D, Lombardero M, Sanz ML. Two different profiles of peach allergy in the north of Spain. *Allergy.* 2007;62:408-14.
8. Hartz C, San Miguel-Moncin Mdel M, Cistero-Bahima A, Fotisch K, Metzner KJ, Fortunato D, Lidholm J, Vieths S, Scheurer S. Molecular characterisation of Lac s 1, the major allergen from lettuce (*Lactuca sativa*). *Mol Immunol.* 2007;44:2820-30.
9. Asero R, Mistrello G, Roncarolo D, Amato S. PT with heat-processed apple peel extract to detect LTP hypersensitivity. *Allerg Immunol (Paris).* 2006;38:351-4.
10. Lauer I, Miguel-Moncin MS, Abel T, Foetisch K, Hartz C, Fortunato D, Cistero-Bahima A, Vieths S, Scheurer S. Identification of a plane pollen lipid transfer protein (Pla a 3) and its immunological relation to the peach lipid-transfer protein, Pru p 3. *Clin Exp Allergy.* 2007;37:261-9.
15. Sanchez-Monge R, Blanco C, Lopez-Torres G, Cumplido J, Recas M, Figueroa J, Carrillo T, Salcedo G. Differential allergen sensitization patterns in chestnut allergy with or without associated latex-fruit syndrome. *J Allergy Clin Immunol.* 2006;118:705-10.
16. Vassilopoulou E, Rigby N, Moreno FJ, Zuidmeer L, Akkerdaas J, Tassios I, Papadopoulos NG, Saxoni-Papageorgiou P, van Ree R, Mills C. Effect of in vitro gastric and duodenal digestion on the allergenicity of grape lipid transfer protein. *J Allergy Clin Immunol.* 2006;118:473-80.
17. Asero R, Mistrello G, Amato S, Roncarolo D, Martinelli A, Zaccarini M. Peach fuzz contains large amounts of lipid transfer protein: is this the cause of the high prevalence of sensitization to LTP in Mediterranean countries? *Allerg Immunol (Paris).* 2006;38:118-21.
18. Reuter A, Lidholm J, Andersson K, Ostling J, Lundberg M, Scheurer S, Enrique E, Cistero-Bahima A, San Miguel-Moncin M, Ballmer-Weber BK, Vieths S. A critical assessment of allergen component-based in vitro diagnosis in cherry allergy across Europe. *Clin Exp Allergy.* 2006;36:815-23.
19. Palacin A, Cumplido J, Figueroa J, Ahrazem O, Sanchez-Monge R, Carrillo T, Salcedo G, Blanco C. Cabbage lipid transfer protein Bra o 3 is a major allergen responsible for cross-reactivity between plant foods and pollens. *J Allergy Clin Immunol.* 2006;117:1423-9.
20. Zuidmeer L, Salentijn E, Rivas MF, Mancebo EG, Asero R, Matos CI, Pelgrom KT, Gilissen LJ, van Ree R. The role of profilin and lipid transfer protein in strawberry allergy in the Mediterranean area. *Clin Exp Allergy.* 2006;36:666-75.
21. Cudowska B, Kaczmarek M, Restani P. Immunoblotting in the diagnosis of cross-reactivity in children allergic to birch. *Rocz Akad Med Bialymst.* 2005;50:268-73.
22. Sancho AI, Rigby NM, Zuidmeer L, Asero R, Mistrello G, Amato S, Gonzalez-Mancebo E, Fernandez-Rivas M, van Ree R, Mills EN. The effect of thermal processing on the IgE reactivity of the non-specific lipid transfer protein from apple, Mal d 3. *Allergy.* 2005;60:1262-8.
23. Ballmer-Weber BK, Wangorsch A, Bohle B, Kaul S, Kundig T, Fotisch K, van Ree R, Vieths S. Component-resolved in vitro diagnosis in carrot allergy: does the use of recombinant carrot allergens improve the reliability of the diagnostic procedure? *Clin Exp Allergy.* 2005;35:970-8.

24. Zuidmeer L, van Leeuwen WA, Budde IK, Cornelissen J, Bulder I, Rafalska I, Besoli NT, Akkerdaas JH, Asero R, Fernandez Rivas M, Gonzalez Mancebo E, van Ree R. Lipid transfer proteins from fruit: cloning, expression and quantification. *Int Arch Allergy Immunol.* 2005;137:273-81.
25. Asero R, Mistrello G, Roncarolo D, Amato S, Falagiani P. Why do lipid transfer protein-hypersensitive patients tolerate bean (and other legumes)? *Int Arch Allergy Immunol.* 2005;137:236-40.
26. Ahrazem O, Ibanez MD, Lopez-Torrejon G, Sanchez-Monge R, Sastre J, Lombardero M, Barber D, Salcedo G. Lipid transfer proteins and allergy to oranges. *Int Arch Allergy Immunol.* 2005;137:201-10.
27. Schad SG, Trcka J, Vieths S, Scheurer S, Conti A, Brocker EB, Trautmann A. Wine anaphylaxis in a German patient: IgE-mediated allergy against a lipid transfer protein of grapes. *Int Arch Allergy Immunol.* 2005;136:159-64.
28. Pastorello EA, Farioli L, Pravettoni V, Robino AM, Scibilia J, Fortunato D, Conti A, Borgonovo L, Bengtsson A, Ortolani C. Lipid transfer protein and vicilin are important walnut allergens in patients not allergic to pollen. *J Allergy Clin Immunol.* 2004;114:908-14.
29. Lombardero M, Garcia-Selles FJ, Polo F, Jimeno L, Chamorro MJ, Garcia-Casado G, Sanchez-Monge R, Diaz-Perales A, Salcedo G, Barber D. Prevalence of sensitization to *Artemisia* allergens Art v 1, Art v 3 and Art v 60 kDa. Cross-reactivity among Art v 3 and other relevant lipid-transfer protein allergens. *Clin Exp Allergy.* 2004;34:1415-21.
30. Akkerdaas JH, Wensing M, Knulst AC, Aalberse RC, Hefle SL, van Ree R. In vitro and in vivo characterization of hazelnut skin prick test extracts. *Arb Paul Ehrlich Inst Bundesamt Sera Impfstoffe Frankf A M.* 2003;(94):87-95; discussion 96.
31. Asero R, Mistrello G, Roncarolo D, Amato S. Relationship between peach lipid transfer protein specific IgE levels and hypersensitivity to non-Rosaceae vegetable foods in patients allergic to lipid transfer protein. *Ann Allergy Asthma Immunol.* 2004;92:268-72.
32. Tabar AI, Alvarez-Puebla MJ, Gomez B, Sanchez-Monge R, Garcia BE, Echechipia S, Olaguibel JM, Salcedo G. Diversity of asparagus allergy: clinical and immunological features. *Clin Exp Allergy.* 2004;34:131-6.
33. Schocker F, Luttkopf D, Scheurer S, Petersen A, Cistero-Bahima A, Enrique E, San Miguel-Moncin M, Akkerdaas J, van Ree R, Vieths S, Becker WM. Recombinant lipid transfer protein Cor a 8 from hazelnut: a new tool for in vitro diagnosis of potentially severe hazelnut allergy. *J Allergy Clin Immunol.* 2004;113:141-7.
34. Fernandez-Rivas M, Gonzalez-Mancebo E, Rodriguez-Perez R, Benito C, Sanchez-Monge R, Salcedo G, Alonso MD, Rosado A, Tejedor MA, Vila C, Casas ML. Clinically relevant peach allergy is related to peach lipid transfer protein, Pru p 3, in the Spanish population. *J Allergy Clin Immunol.* 2003;112:789-95.
35. Pastorello EA, Pompei C, Pravettoni V, Farioli L, Calamari AM, Scibilia J, Robino AM, Conti A, Iametti S, Fortunato D, Bonomi S, Ortolani C. Lipid-transfer protein is the major maize allergen maintaining IgE-binding activity after cooking at 100 degrees C, as demonstrated in anaphylactic patients and patients with positive double-blind, placebo-controlled food challenge results. *J Allergy Clin Immunol.* 2003;112:775-83.
36. Garcia-Casado G, Pacios LF, Diaz-Perales A, Sanchez-Monge R, Lombardero M, Garcia-Selles FJ, Polo F, Barber D, Salcedo G. Identification of IgE-binding epitopes of the major peach allergen Pru p 3. *J Allergy Clin Immunol.* 2003;112:599-605.
37. San Miguel-Moncin M, Krail M, Scheurer S, Enrique E, Alonso R, Conti A, Cistero-Bahima A, Vieths S. Lettuce anaphylaxis: identification of a lipid transfer protein as the major allergen. *Allergy.* 2003 ;58:511-7.
38. Beezhold DH, Hickey VL, Kostyal DA, Puhl H, Zuidmeer L, van Ree R, Sussman GL. Lipid transfer protein from *Hevea brasiliensis* (Hev b 12), a cross-reactive latex protein. *Ann Allergy Asthma Immunol.* 2003;90:439-45.
39. Asturias JA, Gomez-Bayon N, Eserverri JL, Martinez A. Par j 1 and Par j 2, the major allergens from *Parietaria judaica* pollen, have similar immunoglobulin E epitopes. *Clin Exp Allergy.* 2003;33:518-24.
40. Pastorello EA, Farioli L, Pravettoni V, Ortolani C, Fortunato D, Giuffrida MG, Perono Garoffo L, Calamari AM, Brenna O, Conti A. Identification of grape and wine allergens as an endochitinase 4, a lipid-transfer protein, and a thaumatin. *J Allergy Clin Immunol.* 2003;111:350-9.
41. Diaz-Perales A, Tabar AI, Sanchez-Monge R, Garcia BE, Gomez B, Barber D, Salcedo G. Characterization of asparagus allergens: a relevant role of lipid transfer proteins. *J Allergy Clin Immunol.* 2002;110:790-6.
42. Carnes J, Fernandez-Caldas E, Gallego MT, Ferrer A, Cuesta-Herranz J. Pru p 3 (LTP) content in peach extracts. *Allergy.* 2002;57:1071-5.
43. Asero R, Mistrello G, Roncarolo D, Amato S, Caldironi G, Barocci F, van Ree R. Immunological cross-reactivity between lipid transfer proteins from botanically unrelated plant-derived foods: a clinical study. *Allergy.* 2002;57:900-6.
44. Florido Lopez JF, Quiralte Enriquez J, Arias de Saavedra Alias JM, Saenz de San Pedro B, Martin Casanez E. An allergen from *Olea europaea* pollen (Ole e 7) is associated with plant-derived food anaphylaxis. *Allergy.* 2002;57 Suppl 71:53-9.
45. Pastorello EA, Pravettoni V, Farioli L, Rivolta F, Conti A, Ispano M, Fortunato D, Bengtsson A, Bianchi M. Hypersensitivity to mugwort (*Artemisia vulgaris*) in patients with peach allergy is due to a common lipid transfer protein allergen and is often without clinical expression. *J Allergy Clin Immunol.* 2002;110:310-7.
46. Garcia-Selles FJ, Diaz-Perales A, Sanchez-Monge R, Alcantara M, Lombardero M, Barber D, Salcedo G, Fernandez-Rivas M. Patterns of reactivity to lipid transfer proteins of plant foods and *Artemisia* pollen: an in vivo study. *Int Arch Allergy Immunol.* 2002;128:115-22.
47. Pastorello EA, Vieths S, Pravettoni V, Farioli L, Trambaioli C, Fortunato D, Luttkopf D, Calamari M, Ansaloni R, Scibilia J, Ballmer-Weber BK, Poulsen LK, Wutrich B, Hansen KS, Robino AM, Ortolani C, Conti A. Identification of hazelnut major allergens in sensitive patients with positivedouble-blind, placebo-controlled food challenge results. *J Allergy Clin Immunol.* 2002;109:563-70.
48. Asero R, Mistrello G, Roncarolo D, Casarini M, Falagiani P. Allergy to nonspecific lipid transfer proteins in Rosaceae: a comparative study of different in vivo diagnostic methods. *Ann Allergy Asthma Immunol.* 2001;87:68-71.
49. Asero R, Mistrello G, Roncarolo D, Amato S, van Ree R. A case of allergy to beer showing cross-reactivity between lipid transfer proteins. *Ann Allergy Asthma Immunol.* 2001;87:65-7.
50. Pastorello EA, Farioli L, Pravettoni V, Giuffrida MG, Ortolani C, Fortunato D, Trambaioli C, Scibola E, Calamari AM, Robino AM, Conti A. Characterization of the major allergen of plum as a lipid transfer protein. *J Chromatogr B Biomed Sci Appl.* 2001;756:95-103.

51. Conti A, Fortunato D, Ortolani C, Giuffrida MG, Pravettoni V, Napolitano L, Farioli L, Perono Garoffo L, Trambaioli C, Pastorello EA. Determination of the primary structure of two lipid transfer proteins from apricot (*Prunus armeniaca*). *J Chromatogr B Biomed Sci Appl.* 2001;756:123-9.
52. Asero R, Mistrello G, Roncarolo D, de Vries SC, Gautier MF, Ciurana CL, Verbeek E, Mohammadi T, Knul-Brettlova V, Akkerdaas JH, Bulder I, Aalberse RC, van Ree R. Lipid transfer protein: a pan-allergen in plant-derived foods that is highly resistant to pepsin digestion. *Int Arch Allergy Immunol.* 2001;124:67-9.
53. Scheurer S, Pastorello EA, Wangorsch A, Kastner M, Hausteiner D, Vieths S. Recombinant allergens Pru av 1 and Pru av 4 and a newly identified lipid transfer protein in the in vitro diagnosis of cherry allergy. *J Allergy Clin Immunol.* 2001;107:724-31.
54. Pastorello EA, Farioli L, Pravettoni V, Spano M, Scibola E, Trambaioli C, Giuffrida MG, Ansaloni R, Godovac-Zimmermann J, Conti A, Fortunato D, Ortolani C. The maize major allergen, which is responsible for food-induced allergic reactions, is a lipid transfer protein. *J Allergy Clin Immunol.* 2000;106:744-51.
55. Asero R, Mistrello G, Roncarolo D, de Vries SC, Gautier MF, Ciurana CL, Verbeek E, Mohammadi T, Knul-Brettlova V, Akkerdaas JH, Bulder I, Aalberse RC, van Ree R. Lipid transfer protein: a pan-allergen in plant-derived foods that is highly resistant to pepsin digestion. *Int Arch Allergy Immunol.* 2000;122:20-32.
56. Pastorello EA, D'Ambrosio FP, Pravettoni V, Farioli L, Giuffrida G, Monza M, Ansaloni R, Fortunato D, Scibola E, Rivolta F, Incorvaia C, Bengtsson A, Conti A, Ortolani C. Evidence for a lipid transfer protein as the major allergen of apricot. *J Allergy Clin Immunol.* 2000 Feb;105:371-7.
57. Asero R. Detection and clinical characterization of patients with oral allergy syndrome caused by stable allergens in Rosaceae and nuts. *Ann Allergy Asthma Immunol.* 1999;83:377-83.
7. Crespo JF, Retzek M, Foetisch K, Sierra-Maestro E, Cid-Sanchez AB, Pascual CY, Conti A, Feliu A, Rodriguez J, Vieths S, Scheurer S. Germin-like protein Cit s 1 and profilin Cit s 2 are major allergens in orange (*Citrus sinensis*) fruits. *Mol Nutr Food Res.* 2006;50:282-90.
8. Lopez-Torres G, Ibanez MD, Ahrazem O, Sanchez-Monge R, Sastre J, Lombardero M, Barber D, Salcedo G. Isolation, cloning and allergenic reactivity of natural profilin Cit s 2, a major orange allergen. *Allergy.* 2005;60:1424-9.
9. Sankian M, Varasteh A, Pazouki N, Mahmoudi M. Sequence homology: a poor predictive value for profilins cross-reactivity. *Clin Mol Allergy.* 2005;3:13.
10. Lopez-Torres G, Crespo JF, Sanchez-Monge R, Sanchez-Jimenez M, Alvarez J, Rodriguez J, Salcedo G. Allergenic reactivity of the melon profilin Cuc m 2 and its identification as major allergen. *Clin Exp Allergy.* 2005;35:1065-72.
11. Asturias JA, Ibarrola I, Fernandez J, Arilla MC, Gonzalez-Rioja R, Martinez A. Pho d 2, a major allergen from date palm pollen, is a profilin: cloning, sequencing, and immunoglobulin E cross-reactivity with other profilins. *Clin Exp Allergy.* 2005;35:374-81.
12. Enrique E, Alonso R, Bartolome B, San Miguel-Moncin M, Bartra J, Fernandez-Parra B, Tella R, Asturias JA, Ibarrola I, Martinez A, Cistero-Bahima A. IgE reactivity to profilin in *Platanus acerifolia* pollen-sensitized subjects with plant-derived food allergy. *J Investig Allergol Clin Immunol.* 2004;14:335-42.
13. Ghunaim N, Gronlund H, Kronqvist M, Gronneberg R, Soderstrom L, Ahlstedt S, van Hage-Hamsten M. Antibody profiles and self-reported symptoms to pollen-related food allergens in grass pollen-allergic patients from northern Europe. *Allergy.* 2005;60:185-91.
14. Cudowska B, Kaczmarek M. Diagnostic value of birch recombinant allergens (rBet v 1, profilin rBet v 2) in children with pollen-related food allergy. *Rocz Akad Med Bialymst.* 2004;49:111-5.
15. Wagner S, Radauer C, Hafner C, Fuchs H, Jensen-Jarolim E, Wuthrich B, Scheiner O, Breiteneder H. Characterization of cross-reactive bell pepper allergens involved in the latex-fruit syndrome. *Clin Exp Allergy.* 2004;34:1739-46.
16. Scheurer S, Lauer I, Foetisch K, San Miguel Moncin M, Retzek M, Hartz C, Enrique E, Lidholm J, Cistero-Bahima A, Vieths S. Strong allergenicity of Pru av 3, the lipid transfer protein from cherry, is related to high stability against thermal processing and digestion. *J Allergy Clin Immunol.* 2004;114:900-7.
17. Barderas R, Villalba M, Rodriguez R. Recombinant expression, purification and cross-reactivity of chenopod profilin: rChe a 2 as a good marker for profilin sensitization. *Biol Chem.* 2004;385:731-7.
18. Barderas R, Villalba M, Pascual CY, Batanero E, Rodriguez R. Profilin (Che a 2) and polcalcin (Che a 3) are relevant allergens of *Chenopodium album* pollen: isolation, amino acid sequences, and immunologic properties. *J Allergy Clin Immunol.* 2004;113:1192-8.
19. Westphal S, Kempf W, Foetisch K, Retzek M, Vieths S, Scheurer S. Tomato profilin Lyc e 1: IgE cross-reactivity and allergenic potency. *Allergy.* 2004;59:526-32.
20. Asero R, Mistrello G, Roncarolo D, Amato S. Parietaria profilin shows only limited cross-reactivity with birch and grass profilins. *Int Arch Allergy Immunol.* 2004;133:121-4.
21. Ebo DG, Hagendorens MM, Bridts CH, De Clerck LS, Stevens WJ. Sensitization to cross-reactive carbohydrate determinants and the ubiquitous protein profilin: mimickers of allergy. *Clin Exp Allergy.* 2004;34:137-44.

## Literatur zu Profilinen

1. Lopez-Torres G, Diaz-Perales A, Rodriguez J, Sanchez-Monge R, Crespo JF, Salcedo G, Pacios LF. An experimental and modeling-based approach to locate IgE epitopes of plant profilin allergens. *J Allergy Clin Immunol.* 2007;119:1481-8.
2. Gamboa PM, Caceres O, Antepara I, Sanchez-Monge R, Ahrazem O, Salcedo G, Barber D, Lombardero M, Sanz ML. Two different profiles of peach allergy in the north of Spain. *Allergy.* 2007;62:408-14.
3. Ma Y, Zuidmeer L, Bohle B, Bolhaar ST, Gadermaier G, Gonzalez-Mancebo E, Fernandez-Rivas M, Knulst AC, Himly M, Asero R, Ebner C, van Ree R, Ferreira F, Breiteneder H, Hoffmann-Sommergruber K. Characterization of recombinant Mal d 4 and its application for component-resolved diagnosis of apple allergy. *Clin Exp Allergy.* 2007;36:1087-96.
4. Radauer C, Willeroider M, Fuchs H, Hoffmann-Sommergruber K, Thalhammer J, Ferreira F, Scheiner O, Breiteneder H. Cross-reactive and species-specific immunoglobulin E epitopes of plant profilins: an experimental and structure-based analysis. *Clin Exp Allergy.* 2006;36:920-9.
5. Compes E, Hernandez E, Quirce S, Palomares O, Rodriguez R, Cuesta J, Sastre J, Villalba M. Hypersensitivity to black locust (*Robinia pseudoacacia*) pollen: "allergy mirages". *Ann Allergy Asthma Immunol.* 2006;96:586-92.
6. Zuidmeer L, Salentijn E, Rivas MF, Mancebo EG, Asero R, Matos CI, Pelgrom KT, Gilissen LJ, van Ree R. The role of profilin and lipid transfer protein in strawberry allergy in the Mediterranean area. *Clin Exp Allergy.* 2006;36:666-75.

22. De Amici M, Mosca M, Vignini M, Quaglini S, Moratti R. Recombinant birch allergens (Bet v 1 and Bet v 2) and the oral allergy syndrome in patients allergic to birch pollen. *Ann Allergy Asthma Immunol.* 2003;91:490-2.
23. Fernandez-Rivas M, Gonzalez-Mancebo E, Rodriguez-Perez R, Benito C, Sanchez-Monge R, Salcedo G, Alonso MD, Rosado A, Tejedor MA, Vila C, Casas ML. Clinically relevant peach allergy is related to peach lipid transfer protein, Prup 3, in the Spanish population. *J Allergy Clin Immunol.* 2003;112:789-95.
24. Mari A, Wallner M, Ferreira F. Fagales pollen sensitization in a birch-free area: a respiratory cohort survey using Fagales pollen extracts and birch recombinant allergens (rBet v 1, rBet v 2, rBet v 4). *Clin Exp Allergy.* 2003;33:1419-28.
25. Willeroider M, Fuchs H, Ballmer-Weber BK, Focke M, Susani M, Thalhamer J, Ferreira F, Wuthrich B, Scheiner O, Breiteneder H, Hoffmann-Sommergruber K. Cloning and molecular and immunological characterisation of two new food allergens, Cap a 2 and Lyc e 1, profilins from bell pepper (*Capsicum annuum*) and Tomato (*Lycopersicon esculentum*). *Int Arch Allergy Immunol.* 2003;131:245-55.
26. Rossi RE, Monasterolo G, Monasterolo S. Detection of specific IgE antibodies in the sera of patients allergic to birch pollen using recombinant allergens Bet v 1, Bet v 2, Bet v 4: evaluation of different IgE reactivity profiles. *Allergy.* 2003;58:929-32.
27. Asero R, Mistrello G, Roncarolo D, Amato S, Zanoni D, Barocci F, Caldironi G. Detection of clinical markers of sensitization to profilin in patients allergic to plant-derived foods. *J Allergy Clin Immunol.* 2003;112:427-32.
28. Rodriguez-Perez R, Fernandez-Rivas M, Gonzalez-Mancebo E, Sanchez-Monge R, Diaz-Perales A, Salcedo G. Peach profilin: cloning, heterologous expression and cross-reactivity with Bet v 2. *Allergy.* 2003;58:635-40.
29. Rodriguez-Perez R, Crespo JF, Rodriguez J, Salcedo G. Profilin is a relevant melon allergen susceptible to pepsin digestion in patients with oral allergy syndrome. *J Allergy Clin Immunol.* 2003;111:634-9.
30. Diaz-Perales A, Tabar AI, Sanchez-Monge R, Garcia BE, Gomez B, Barber D, Salcedo G. Characterization of asparagus allergens: a relevant role of lipid transfer proteins. *J Allergy Clin Immunol.* 2002;110:790-6.
31. Wensing M, Akkerdaas JH, van Leeuwen WA, Stapel SO, Bruijnzeel-Koomen CA, Aalberse RC, Bast BJ, Knulst AC, van Ree R. IgE to Bet v 1 and profilin: cross-reactivity patterns and clinical relevance. *J Allergy Clin Immunol.* 2002;110:435-42.
32. Florido Lopez JF, Quiralte Enriquez J, Arias de Saavedra Alias JM, Saenz de San Pedro B, Martin Casanez E. An allergen from *Olea europaea* pollen (Ole e 7) is associated with plant-derived food anaphylaxis. *Allergy.* 2002;57 Suppl 71:53-9.
33. Nieto A, Mazon A, Boquete M, Carballada F, Asturias JA, Martinez J, Martinez A. Assessment of profilin as an allergen for latex-sensitized patients. *Allergy.* 2002;57:776-84.
34. Ballmer-Weber BK, Scheurer S, Fritsche P, Enrique E, Cistero-Bahima A, Haase T, Wuthrich B. Component-resolved diagnosis with recombinant allergens in patients with cherry allergy. *J Allergy Clin Immunol.* 2002;110:167-73.
35. Reindl J, Rihs HP, Scheurer S, Wangorsch A, Haustein D, Vieths S. IgE reactivity to profilin in pollen-sensitized subjects with adverse reactions to banana and pineapple. *Int Arch Allergy Immunol.* 2002;128:105-14.
36. Niederberger V, Purohit A, Oster JP, Spitzauer S, Valenta R, Pauli G. The allergen profile of ash (*Fraxinus excelsior*) pollen: cross-reactivity with allergens from various plant species. *Clin Exp Allergy.* 2002;32:933-41.
37. Kwaasi AA, Harfi HA, Parhar RS, Saleh S, Collison KS, Panzani RC, Al-Sedairy ST, Al-Mohanna FA. Cross-reactivities between date palm (*Phoenix dactylifera* L.) polypeptides and foods implicated in the oral allergy syndrome. *Allergy.* 2002;57:508-18.
38. Ballmer-Weber BK, Hoffmann A, Wuthrich B, Luttkopf D, Pompei C, Wangorsch A, Kastner M, Vieths S. Influence of food processing on the allergenicity of celery: DBPCFC with celery spice and cooked celery in patients with celery allergy. *Allergy.* 2002;57:228-35.
39. Ballmer-Weber BK, Wuthrich B, Wangorsch A, Fotisch K, Altmann F, Vieths S. Carrot allergy: double-blinded, placebo-controlled food challenge and identification of allergens. *J Allergy Clin Immunol.* 2001;108:301-7.
40. Ganglberger E, Radauer C, Wagner S, Riordain G, Beezhold DH, Brehler R, Niggemann B, Scheiner O, Jensen-Jarolim E, Breiteneder H. Hev b 8, the *Hevea brasiliensis* latex profilin, is a cross-reactive allergen of latex, plant foods and pollen. *Int Arch Allergy Immunol.* 2001;125:216-27.
41. Benitez D, Garcia-Ortega P, Picado C, Mila J, Vives J, Martinez J, Vilella R. Specific immune response to *Phleum pratense* plant profilin in atopic patients and control subjects. *Allergol Immunopathol (Madr).* 2001;29:9-15.
42. Scheurer S, Wangorsch A, Nerkamp J, Skov PS, Ballmer-Weber B, Wuthrich B, Haustein D, Vieths S. Cross-reactivity within the profilin panallergen family investigated by comparison of recombinant profilins from pear (Pyr c 4), cherry (Pru av 4) and celery (Api g 4) with birch pollen profilin Bet v 2. *J Chromatogr B Biomed Sci Appl.* 2001;756:315-25.
43. Mari A. Multiple pollen sensitization: a molecular approach to the diagnosis. *Int Arch Allergy Immunol.* 2001;125:57-65.
44. Asero R, Mistrello G, Roncarolo D, de Vries SC, Gautier MF, Ciurana CL, Verbeek E, Mohammadi T, Knul-Brettlova V, Akkerdaas JH, Bulder I, Aalberse RC, van Ree R. Lipid transfer protein: a pan-allergen in plant-derived foods that is highly resistant to pepsin digestion. *Int Arch Allergy Immunol.* 2001;124:67-9.
45. Anliker MD, Reindl J, Vieths S, Wuthrich B. Allergy caused by ingestion of persimmon (*Diospyros kaki*): detection of specific IgE and cross-reactivity to profilin and carbohydrate determinants. *J Allergy Clin Immunol.* 2001;107:718-23.
46. Hemmer W, Focke M, Wantke F, Gotz M, Jarisch R, Jager S, Gotz M. Ash (*Fraxinus excelsior*)-pollen allergy in central Europe: specific role of pollen panallergens and the major allergen of ash pollen, Fra e 1. *Allergy.* 2000;55:923-30.
47. Luttkopf D, Ballmer-Weber BK, Wuthrich B, Vieths S. Celery allergens in patients with positive double-blind placebo-controlled food challenge. *J Allergy Clin Immunol.* 2000;106:390-9.
48. Scheurer S, Wangorsch A, Haustein D, Vieths S. Cloning of the minor allergen Api g 4 profilin from celery (*Apium graveolens*) and its cross-reactivity with birch pollen profilin Bet v 2. *Clin Exp Allergy.* 2000;30:962-71.
49. Ganglberger E, Radauer C, Grimm R, Hoffmann-Sommergruber K, Breiteneder H, Scheiner O, Jensen-Jarolim E. N-terminal sequences of high molecular weight allergens from celery tuber. *Clin Exp Allergy.* 2000;30:566-70.
50. Diez-Gomez ML, Quirce S, Cuevas M, Sanchez-Fernandez C, Baz G, Moradiellos FJ, Martinez A. Fruit-pollen-latex cross-reactivity: implication of profilin (Bet v 2). *Allergy.* 1999;54:951-61.

51. Jensen-Jarolim E, Gerstmayer G, Kraft D, Scheiner O, Ebner H, Ebner C. Serological characterization of allergens in poppy seeds. *Clin Exp Allergy*. 1999;29:1075-9.
52. Martinez A, Asturias JA, Palacios R, Sanz ML, Sanchez G, Oehling A, Martinez J. Identification of a 36-kDa olive-pollen allergen by in vitro and in vivo studies. *Allergy*. 1999;54:584-92.
53. Van Ree R, Van Leeuwen WA, Akkerdaas JH, Aalberse RC. How far can we simplify in vitro diagnostics for Fagales tree pollen allergy? A study with three whole pollen extracts and purified natural and recombinant allergens. *Clin Exp Allergy*. 1999;29:848-55.
54. Niederberger V, Pauli G, Gronlund H, Froschl R, Rumpold H, Kraft D, Valenta R, Spitzauer S. Recombinant birch pollen allergens (rBet v 1 and rBet v 2) contain most of the IgE epitopes present in birch, alder, hornbeam, hazel, and oak pollen: quantitative IgE inhibition study with sera from different populations. *J Allergy Clin Immunol*. 1998;102:579-91.
55. Asturias JA, Arilla MC, Gomez-Bayon N, Aguirre M, Martinez A, Palacios R, Martinez J. Cloning and immunological characterization of the allergen Hel a 2 (profilin) from sunflower pollen. *Mol Immunol*. 1998;35:469-78.
56. Focke M, Hemmer W, Hayek B, Gotz M, Jarisch R. Identification of allergens in oilseed rape (*Brassica napus*) pollen. *Int Arch Allergy Immunol*. 1998;117:105-12.
57. Daschner A, Crespo JF, Pascual CY. Specific IgE to recombinant vegetal panallergen (rBet v 2) and fruit allergy in pollinic patients. *Allergy*. 1998;53:614-8.
58. Jensen-Jarolim E, Santner B, Leitner A, Grimm R, Scheiner O, Ebner C, Breiteneder H. Bell peppers (*Capsicum annuum*) express allergens (profilin, pathogenesis-related protein P23 and Bet v 1) depending on the horticultural strain. *Int Arch Allergy Immunol*. 1998;116:103-9.
59. Rudeschko O, Fahlbusch B, Steurich F, Schlenvoigt G, Jager L. Kiwi allergens and their cross-reactivity with birch, rye, timothy, and mugwort pollen. *J Investig Allergol Clin Immunol*. 1998;8:78-84.
60. Vallverdu A, Asturias JA, Arilla MC, Gomez-Bayon N, Martinez A, Martinez J, Palacios R. Characterization of recombinant *Mercurialis annua* major allergen Mer a 1 (profilin). *J Allergy Clin Immunol*. 1998;101:363-70.
61. Niederberger V, Laffer S, Froschl R, Kraft D, Rumpold H, Kapiotis S, Valenta R, Spitzauer S. IgE antibodies to recombinant pollen allergens (Phl p 1, Phl p 2, Phl p 5, and Bet v 2) account for a high percentage of grass pollen-specific IgE. *J Allergy Clin Immunol*. 1998;101:258-64.
62. Hirschehr R, Heppner C, Spitzauer S, Sperr WR, Valent P, Berger U, Horak F, Jager S, Kraft D, Valenta R. Identification of common allergenic structures in mugwort and ragweed pollen. *J Allergy Clin Immunol*. 1998;10:196-206.
63. Leitner A, Jensen-Jarolim E, Grimm R, Wuthrich B, Ebner H, Scheiner O, Kraft D, Ebner C. Allergens in pepper and paprika. Immunologic investigation of the celery-birch-mugwort-spice syndrome. *Allergy*. 1998;53:36-41.
64. Fernandez-Rivas M, van Ree R, Cuevas M. Allergy to Rosaceae fruits without related pollinosis. *J Allergy Clin Immunol*. 1997;100:728-33.
65. Jensen-Jarolim E, Leitner A, Hirschehr R, Kraft D, Wuthrich B, Scheiner O, Graf J, Ebner C. Characterization of allergens in Apiaceae spices: anise, fennel, coriander and cumin. *Clin Exp Allergy*. 1997;27:1299-306.
66. Asturias JA, Arilla MC, Gomez-Bayon N, Martinez J, Martinez A, Palacios R. Cloning and expression of the panallergen profilin and the major allergen (Ole e 1) from olive tree pollen. *J Allergy Clin Immunol*. 1997;100:365-72.
67. Fuchs T, Spitzauer S, Vente C, Hevler J, Kapiotis S, Rumpold H, Kraft D, Valenta R. Natural latex, grass pollen, and weed pollen share IgE epitopes. *J Allergy Clin Immunol*. 1997;100:356-64.
68. Elfman L, Svensson M, Lidholm J, Pauli G, Valenta R. Different profiles in specific IgE to rBet v 1 and rBet v 2 in patients allergic to birch pollen from six countries. *Int Arch Allergy Immunol*. 1997;113:249-51.
69. Rossi RE, Monasterolo G, Operti D, Corsi M. Evaluation of recombinant allergens Bet v 1 and Bet v 2 (profilin) by Pharmacia CAP system in patients with pollen-related allergy to birch and apple. *Allergy*. 1996;51:940-5.
70. Bauer L, Ebner C, Hirschehr R, Wuthrich B, Pichler C, Fritsch R, Scheiner O, Kraft D. IgE cross-reactivity between birch pollen, mugwort pollen and celery is due to at least three distinct cross-reacting allergens: immunoblot investigation of the birch-mugwort-celery syndrome. *Clin Exp Allergy*. 1996;26:1161-70.
71. Laffer S, Spitzauer S, Susani M, Pairleitner H, Schweiger C, Gronlund H, Menz G, Pauli G, Ishii T, Nolte H, Ebner C, Sehon AH, Kraft D, Eichler HG, Valenta R. Comparison of recombinant timothy grass pollen allergens with natural extract for diagnosis of grass pollen allergy in different populations. *J Allergy Clin Immunol*. 1996;98:652-8.
72. Pauli G, Oster JP, Deviller P, Heiss S, Bessot JC, Susani M, Ferreira F, Kraft D, Valenta R. Skin testing with recombinant allergens rBet v 1 and birch profilin, rBet v 2: diagnostic value for birch pollen and associated allergies. *J Allergy Clin Immunol*. 1996;97:1100-9.
73. Eng PA, Yman L, Maaninen E, Wuthrich B. Inhalant allergy to fresh asparagus. *Clin Exp Allergy*. 1996;26:330-4.
74. Fah J, Wuthrich B, Vieths S. Anaphylactic reaction to lychee fruit: evidence for sensitization to profilin. *Clin Exp Allergy*. 1995;25:1018-23.
75. Ebner C, Hirschehr R, Bauer L, Breiteneder H, Valenta R, Ebner H, Kraft D, Scheiner O. Identification of allergens in fruits and vegetables: IgE cross-reactivities with the important birch pollen allergens Bet v 1 and Bet v 2 (birch profilin). *J Allergy Clin Immunol*. 1995;95:962-9.
76. Vallier P, Baland S, Harf R, Valenta R, Deviller P. Identification of profilin as an IgE-binding component in latex from *Hevea brasiliensis*: clinical implications. *Clin Exp Allergy*. 1995;25:332-9.
77. van Ree R, Fernandez-Rivas M, Cuevas M, van Wijngaarden M, Aalberse RC. Pollen-related allergy to peach and apple: an important role for profilin. *J Allergy Clin Immunol*. 1995;95:726-34.
78. Valenta R, Duchene M, Ebner C, Valent P, Sillaber C, Deviller P, Ferreira F, Tejkl M, Edelmann H, Kraft D, et al. Profilins constitute a novel family of functional plant pan-allergens. *J Exp Med*. 1992;175:377-85.
79. Valenta R, Vrtala S, Ebner C, Kraft D, Scheiner O. Diagnosis of grass pollen allergy with recombinant timothy grass (*Phleum pratense*) pollen allergens. *Int Arch Allergy Immunol*. 1992;97:287-94.

## Literatur zu Speicherproteinen

1. Choi SY, Sohn JH, Lee YW, Lee EK, Hong CS, Park JW. Application of the 16-kDa buckwheat 2 S storage albumin protein for diagnosis of clinical reactivity. *Ann Allergy Asthma Immunol.* 2007;99:254-60.
2. Flinterman AE, van Hoffen E, den Hartog Jager CF, Koppelman S, Pasmans SG, Hoekstra MO, Bruijnzeel-Koomen CA, Knulst AC, Knol EF. Children with peanut allergy recognize predominantly Ara h2 and Ara h6, which remains stable over time. *Clin Exp Allergy.* 2007;37:1221-8.
3. McDermott RA, Porterfield HS, El Mezayen R, Burks AW, Pons L, Schlichting DG, Solomon B, Redzic JS, Harbeck RJ, Duncan MW, Hansen KC, Dreskin SC. Contribution of Ara h 2 to peanut-specific, immunoglobulin E-mediated, cell activation. *Clin Exp Allergy.* 2007;37:752-63.
4. Nolan RC, Richmond P, Prescott SL, Mallon DF, Gong G, Franzmann AM, Naidoo R, Loh RK. Skin prick testing predicts peanut challenge outcome in previously allergic or sensitized children with low serum peanut-specific IgE antibody concentration. *Pediatr Allergy Immunol.* 2007;18:224-30.
5. Wainstein BK, Yee A, Jelley D, Ziegler M, Ziegler JB. Combining skin prick, immediate skin application and specific-IgE testing in the diagnosis of peanut allergy in children. *Pediatr Allergy Immunol.* 2007;18:231-9.
6. de Leon MP, Drew AC, Glaspole IN, Suphioglu C, O'Hehir RE, Rolland JM. IgE cross-reactivity between the major peanut allergen Ara h 2 and tree nut allergens. *Mol Immunol.* 2007;44:463-71.
7. Peeters KA, Koppelman SJ, van Hoffen E, van der Tas CW, den Hartog Jager CF, Penninks AH, Hefle SL, Bruijnzeel-Koomen CA, Knol EF, Knulst AC. Does skin prick test reactivity to purified allergens correlate with clinical severity of peanut allergy? *Clin Exp Allergy.* 2007;37:108-15.
8. Andersson K, Ballmer-Weber BK, Cistero-Bahima A, Ostling J, Lauer I, Vieths S, Lidholm J. Enhancement of hazelnut extract for IgE testing by recombinant allergen spiking. *Allergy.* 2007;62:897-904.
9. Benito C, Gonzalez-Mancebo E, de Durana MD, Tolon RM, Fernandez-Rivas M. Identification of a 7S globulin as a novel coconut allergen. *Ann Allergy Asthma Immunol.* 2007;98:580-4.
10. Wallowitz ML, Chen RJ, Tzen JT, Teuber SS. Ses i 6, the sesame 11S globulin, can activate basophils and shows cross-reactivity with walnut in vitro. *Clin Exp Allergy.* 2007;37:929-38.
11. Astier C, Morisset M, Roitel O, Codreanu F, Jacquenet S, Franck P, Ogier V, Petit N, Proust B, Moneret-Vautrin DA, Burks AW, Bihain B, Sampson HA, Kanny G. Predictive value of skin prick tests using recombinant allergens for diagnosis of peanut allergy. *J Allergy Clin Immunol.* 2006;118:250-6.
12. Lin J, Shewry PR, Archer DB, Beyer K, Niggemann B, Haas H, Wilson P, Alcocer MJ. The potential allergenicity of two 2S albumins from soybean (*Glycine max*): a protein microarray approach. *Int Arch Allergy Immunol.* 2006;141:91-102.
13. Leduc V, Moneret-Vautrin DA, Tzen JT, Morisset M, Guerin L, Kanny G. Identification of oleosins as major allergens in sesame seed allergic patients. *Allergy.* 2006;61:349-56.
14. de Leon MP, Drew AC, Glaspole IN, Suphioglu C, Rolland JM, O'Hehir RE. Functional analysis of cross-reactive immunoglobulin E antibodies: peanut-specific immunoglobulin E sensitizes basophils to tree nut allergens. *Clin Exp Allergy.* 2005;35:1056-64.
15. Roberts G, Lack G. Diagnosing peanut allergy with skin prick and specific IgE testing. *J Allergy Clin Immunol.* 2005;115:1291-6.
16. Palomares O, Cuesta-Herranz J, Vereda A, Sirvent S, Villalba M, Rodriguez R. Isolation and identification of an 11S globulin as a new major allergen in mustard seeds. *Ann Allergy Asthma Immunol.* 2005;94:586-92.
17. Robotham JM, Wang F, Seamon V, Teuber SS, Sathe SK, Sampson HA, Beyer K, Seavy M, Roux KH. Ana o 3, an important cashew nut (*Anacardium occidentale* L.) allergen of the 2S albumin family. *J Allergy Clin Immunol.* 2005;115:1284-90.
18. Palomares O, Cuesta-Herranz J, Rodriguez R, Villalba M. A recombinant precursor of the mustard allergen Sin a 1 retains the biochemical and immunological features of the heterodimeric native protein. *Int Arch Allergy Immunol.* 2005;137:18-26.
19. Mittag D, Akkerdaas J, Ballmer-Weber BK, Vogel L, Wensing M, Becker WM, Koppelman SJ, Knulst AC, Helbling A, Hefle SL, Van Ree R, Vieths S. Ara h 8, a Bet v 1-homologous allergen from peanut, is a major allergen in patients with combined birch pollen and peanut allergy. *J Allergy Clin Immunol.* 2004;114:1410-7.
20. Koppelman SJ, Wensing M, Ertmann M, Knulst AC, Knol EF. Relevance of Ara h1, Ara h2 and Ara h3 in peanut-allergic patients, as determined by immunoglobulin E Western blotting, basophil-histamine release and intracutaneous testing: Ara h2 is the most important peanut allergen. *Clin Exp Allergy.* 2004;34:583-90.
21. Comstock SS, McGranahan G, Peterson WR, Teuber SS. Extensive in vitro cross-reactivity to seed storage proteins is present among walnut (*Juglans*) cultivars and species. *Clin Exp Allergy.* 2004;34:1583-90.
22. Wensing M, Knulst AC, Piersma S, O'Kane F, Knol EF, Koppelman SJ. Patients with anaphylaxis to pea can have peanut allergy caused by cross-reactive IgE to vicilin (Ara h 1). *J Allergy Clin Immunol.* 2003;111:420-4.
23. Lopez-Torres G, Salcedo G, Martin-Esteban M, Diaz-Perales A, Pascual CY, Sanchez-Monge R. Len c 1, a major allergen and vicilin from lentil seeds: protein isolation and cDNA cloning. *J Allergy Clin Immunol.* 2003;112:1208-15.
24. Wang F, Robotham JM, Teuber SS, Sathe SK, Roux KH. Ana o 2, a major cashew (*Anacardium occidentale* L.) nut allergen of the legumin family. *Int Arch Allergy Immunol.* 2003;132:27-39.
25. Rancé F, Abbal M, Lauwers-Cancès V. Improved screening for peanut allergy by the combined use of skin prick tests and specific IgE assays. *J Allergy Clin Immunol.* 2002;109:1027-30.
26. Teuber SS, Sathe SK, Peterson WR, Roux KH. Characterization of the soluble allergenic proteins of cashew nut (*Anacardium occidentale* L.). *J Agric Food Chem.* 2002;50:6543-9.
27. Beyer K, Grishina G, Bardina L, Grishin A, Sampson HA. Identification of an 11S globulin as a major hazelnut food allergen in hazelnut-induced systemic reactions. *J Allergy Clin Immunol.* 2002;110:517-23.
28. Wang F, Robotham JM, Teuber SS, Tawde P, Sathe SK, Roux KH. Ana o 1, a cashew (*Anacardium occidentale*) allergen of the vicilin seed storage protein family. *J Allergy Clin Immunol.* 2002;110:160-6.
29. Beyer K, Bardina L, Grishina G, Sampson HA. Identification of sesame seed allergens by 2-dimensional proteomics and Edman sequencing: seed storage proteins as common food allergens. *J Allergy Clin Immunol.* 2002;110:154-9.
30. Robotham JM, Teuber SS, Sathe SK, Roux KH. Linear IgE epitope mapping of the English walnut (*Juglans regia*) major food allergen, Jug r 1. *J Allergy Clin Immunol.* 2002;109:143-9.

31. Roux KH, Teuber SS, Robotham JM, Sathe SK. Detection and stability of the major almond allergen in foods. *J Agric Food Chem.* 2001;49:2131-6.
32. Beardslee TA, Zeece MG, Sarath G, Markwell JP. Soybean glycinin G1 acidic chain shares IgE epitopes with peanut allergen Ara h 3. *Int Arch Allergy Immunol.* 2000;123:299-307.
33. Helm RM, Cockrell G, Connaughton C, Sampson HA, Bannon GA, Beilinson V, Livingstone D, Nielsen NC, Burks AW. A soybean G2 glycinin allergen. 1. Identification and characterization. *Int Arch Allergy Immunol.* 2000;123:205-12.
34. Kelly JD, Hefle SL. 2S methionine-rich protein (SSA) from sunflower seed is an IgE-binding protein. *Allergy.* 2000;55(6):556-60.
35. Seppala U, Alenius H, Turjanmaa K, Reunala T, Palosuo T, Kalkkinen N. Identification of patatin as a novel allergen for children with positive skin prick test responses to raw potato. *J Allergy Clin Immunol.* 1999;103:165-71.
36. Monsalve RI, Gonzalez de la Pena MA, Lopez-Otin C, Fiandor A, Fernandez C, Villalba M, Rodriguez R. Detection, isolation and complete amino acid sequence of an aeroallergenic protein from rapeseed flour. *Clin Exp Allergy.* 1997;27:833-41.
37. Bartolome B, Mendez JD, Armentia A, Vallverdu A, Palacios R. Allergens from Brazil nut: immunochemical characterization. *Allergol Immunopathol (Madr).* 1997;25:135-44.
38. Varjonen E, Bjorksten F, Savolainen J. Stability of cereal allergens. *Clin Exp Allergy.* 1996;26:436-43.
8. Ebo DG, Hagendorens MM, Bridts CH, De Clerck LS, Stevens WJ. Sensitization to cross-reactive carbohydrate determinants and the ubiquitous protein profilin: mimickers of allergy. *Clin Exp Allergy.* 2004;34:137-44.
9. Foetisch K, Westphal S, Lauer I, Retzek M, Altmann F, Kolarich D, Scheurer S, Vieths S. Biological activity of IgE specific for cross-reactive carbohydrate determinants. *J Allergy Clin Immunol.* 2003;111:889-96.
10. Vieths S, Luttkopf D, Reindl J, Anliker MD, Wuthrich B, Ballmer-Weber BK. Allergens in celery and zucchini. *Allergy.* 2002;57 Suppl 72:100-5.
11. Ballmer-Weber BK, Hoffmann A, Wuthrich B, Luttkopf D, Pompei C, Wangorsch A, Kastner M, Vieths S. Influence of food processing on the allergenicity of celery: DBPCFC with celery spice and cooked celery in patients with celery allergy. *Allergy.* 2002;57:228-35.
12. Luttkopf D, Ballmer-Weber BK, Wuthrich B, Vieths S. Celery allergens in patients with positive double-blind placebo-controlled food challenge. *J Allergy Clin Immunol.* 2000;106:390-9.
13. Fotisch K, Altmann F, Hausteiner D, Vieths S. Involvement of carbohydrate epitopes in the IgE response of celery-allergic patients. *Int Arch Allergy Immunol.* 1999;120:30-42.
14. Fotisch K, Fah J, Wuthrich B, Altmann F, Hausteiner D, Vieths S. IgE antibodies specific for carbohydrates in a patient allergic to gum Arabic (*Acacia senegal*). *Allergy.* 1998;53:1043-51.
15. van der Veen MJ, van Ree R, Aalberse RC, Akkerdaas J, Koppelman SJ, Jansen HM, van der Zee JS. Poor biologic activity of cross-reactive IgE directed to carbohydrate determinants of glycoproteins. *J Allergy Clin Immunol.* 1997;100:327-34.

## Literatur zu CCDs

1. Jin C, Hantusch B, Hemmer W, Stadlmann J, Altmann F. Affinity of IgE and IgG against cross-reactive carbohydrate determinants on plant and insect glycoproteins. *J Allergy Clin Immunol.* 2008;121:185-90.
2. Malandain H, Giroux F, Cano Y. The influence of carbohydrate structures present in common allergen sources on specific IgE results. *Eur Ann Allergy Clin Immunol.* 2007;39:216-20.
3. Mahler V, Gutgesell C, Valenta R, Fuchs T. Natural rubber latex and hymenoptera venoms share immunoglobulin E-epitopes accounting for cross-reactive carbohydrate determinants. *Clin Exp Allergy.* 2006;36:1446-56.
4. Jappe U, Raulf-Heimsoth M, Hoffmann M, Burow G, Hubsch-Muller C, Enk A. In vitro hymenoptera venom allergy diagnosis: improved by screening for cross-reactive carbohydrate determinants and reciprocal inhibition. *Allergy.* 2006;61:1220-9.
5. Ballmer-Weber BK, Wangorsch A, Bohle B, Kaul S, Kundig T, Fotisch K, van Ree R, Vieths S. Component-resolved in vitro diagnosis in carrot allergy: does the use of recombinant carrot allergens improve the reliability of the diagnostic procedure? *Clin Exp Allergy.* 2005;35:970-8.
6. Kochuyt AM, Van Hoeyveld EM, Stevens EA. Prevalence and clinical relevance of specific immunoglobulin E to pollen caused by sting-induced specific immunoglobulin E to cross-reacting carbohydrate determinants in Hymenoptera venoms. *Clin Exp Allergy.* 2005;35:441-7.
7. Malandain H. Widening sensitization spectrum through carbohydrate panepitopes – a hypothesis. *Allerg Immunol (Paris).* 2004;36:297-9.

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© Phadia AB, 2008
Rechtmäßiger Hersteller: Phadia AB, Uppsala, Schweden
**Phadia GmbH**, Munzinger Str. 7, D-79111 Freiburg, Tel. +49 761 47 805 0, Fax +49 761 47 805 338
**Phadia Austria GmbH**, Donau-City-Str. 1, A-1220 Wien, Tel. +43 1 270 20 20, Fax +43 1 270 20 20 20
**Phadia AG**, Sennweidstr. 46, CH-6312 Steinhausen, Tel. +41 43 343 40 50, Fax +41 43 343 40 51