HÜHNEREI

Molekulare Allergiediagnostik

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Literatur: 1. Ando H, Moverare R, Kondo Y, Tsuge I, Tanaka A, Borres MP, Urisu A. Utility of ovomucoid-specific IgE concentrations in predicting symptomatic egg allergy. *J Allergy Clin Immunol* 2008;122:583–8. 2. Lemon-Mule H, Sampson HA, Sicherer SH, Shreffler WG, Noone S, Nowak-Węgrzyn A. Immunologic changes in children with egg allergy ingesting extensively heated egg. *J Allergy Clin Immunol* 2008;122:977–83. 3. Urisu A, Ando H, Morita Y, Wada E, Yasaki T, Yamada K, Komada K, Torii S, Goto M, Wakamatsu T. Allergenic activity of heated and ovomucoid-depleted egg white. *J Allergy Clin Immunol* 1997;100:171–6. 4. Bernhisel-Broadbent J, Dintzis HM, Dintzis RZ, Sampson HA. Allergenicity and antigenicity of chicken egg ovomucoid (Gal d III) compared with ovalbumin (Gal d I) in children with egg allergy and in mice. *J Allergy Clin Immunol* 1994;93:1047–59. 5. Montesinos E, Martorell A, Félix R, Cerdá JC. Egg white specific IgE levels in serum as clinical reactivity predictors in the course of egg allergy follow up. *Pediatr Allergy Immunol* 2010;21:634–639. 6. Clarke AT, Skypala I, Leech SC, Ewan PW, Dugué P, Brathwaite N, Huber PAJ, Nasser SM. British Society for Allergy and Clinical Immunology guidelines for the management of egg allergy. *Clin Exp Allergy* 2010;40:1116–1129. 7. Frémont S, Kanny G, Nicolas JP, Moneret-Vautrin DA. Prevalence of lysozyme sensitization in an egg-allergic population. *Allergy* 1997;2:224–228. 8. Pérez-Calderón R, Gonzalo-Garjón MA, Lamilla-Yerga A, Mangas-Santos R, Moreno-Gastón I. Recurrent Angioedema Due to Lysozyme Allergy. *J Investig Allergol Clin Immunol* 2007;17(4):264–266. 9. Jarvinen KM, Beyer K, Vila L, Bardina L, Mishoe M, Sampson HA. Specificity of IgE antibodies to sequential epitopes of hen's egg ovomucoid as a marker for persistence of egg allergy. *Allergy* 2007;62:758–65. 10. Benhamou AH, Caubet JC, Eigenmann PA, Nowak-Węgrzyn, Marcos CP, Reche M and Urisu A. State of the art and new horizons in the diagnosis and management of egg allergy. *Allergy* 2010;65:283–289. 11. Mine Y, Zhang JW. Comparative Study on Antigenicity and Allergenicity of native and denatured egg white proteins. *J Agric Food Chem* 2002;50:2679–2683. 12. Everberg H, Brostedt P, Oman H, Bohman S, Moverare R. Affinity purification of egg-white allergens for improved component-resolved diagnostics. *Int Arch Allergy Immunol* 2011;154:33–41. 13. Quirce S, Marañón F, Umpiérrez A, de las Heras M, Fernández-Caldas E, Sastre J. Chicken serum albumin (Gal d 5) is a partially heat-labile inhalant and food allergen implicated in the bird-egg syndrome. *Allergy*. 2001;56:754–62.

thermoscientific.com/phadia/de

© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.
Legal Manufacturer: Phadia AB, Uppsala, Sweden

Phadia GmbH, Munzinger Str. 7, 79111 Freiburg / Deutschland, Tel. +49 761 47805-0
Phadia Austria GmbH, Donau-City-Str. 1, 1220 Wien / Österreich, Tel. +43 1 2702020
Phadia AG, Senneweidstr. 46, 6312 Steinhausen / Schweiz, Tel. +41 43 3434050

84210260 07/2012

Thermo
SCIENTIFIC

Part of Thermo Fisher Scientific